



I N T E L L I G E N T P R O S T H E T I C D E N T I S T R Y

ENGLISH

CATALOGUE IPD
2025







ipd was founded and originated in the dental sector, specializing in the manufacturing and marketing of PROSTHETIC ABUTMENTS for DENTAL IMPLANTS. Available in more than thirty countries, **ipd** offers a wide range of superior quality products, providing innovative solutions in the field of dental implants/prosthetics. **ipd** is now a global leader in compatible abutments.



Who we are

About us

Our products are manufactured utilizing advanced technology and adhering to the highest quality standards. This includes constant supervision of the exact geometric dimensions and hygiene requirements for the abutments. The components are manufactured under the strictest parameters and procedures established by our engineers in close collaboration with recognised professionals and dentists.

A little background information

We have been operating in Spain since 2004, offering the widest range of dental abutments compatible with the main implant brands. Our goal is to offer technically advanced products of the highest quality and provide the best solutions to our customers.

Quality assurance

At , the quality and excellence in our products is of highest priority.

We work with rigorous quality controls that allow us to offer functional, safe and quality products to meet the most demanding market expectations.

, its management system and its products comply with the applicable regulatory requirements:



- **Medical Device Manufacturing License**

Issued by the AEMPS (Spanish Agency of Medicinal Products and Medical Devices).



- **CE MARK**

All our products have the CE mark according to Regulation EU 2017/745.

Notified Body: IMQ ISTITUTO ITALIANO DEL MARCHIO DI QUALITÀ S.P.A. NB No. 0051.



- **MDSAP (Medical Device Single Audit Program)**

The Quality Management System complies with the requirements of ISO 13485:2016 and of the requirements of the following regulatory authorities: Australia, Brazil, Canada and United States.



- **USA**

IPD is registered in the FDA's Establishment and Device Listing database, with Registration Number 3016954909.



- **Canada**

IPD has the Medical Device Establishment License from Health Canada, MDEL No. 18296.



- **ISO 13485:2016. Medical devices. Quality management systems. Requirements for regulatory purposes.** Certificate No. 9124. IPDE.



- **ISO 9001:2015. Quality management systems. Requirements** Certificate No. 9120. IPDT.



Be Innovative

Custom Interface Ti-Base System

Flexibility made Ti-Base

The most flexible and versatile Ti-Base on the market

Choosing a suitable Ti-Base for each clinical case and restoration is essential for a long term success

Custom Interface Ti-Base system Get to know its innovative design

- + Asymmetrical **tri-lobe engaging** design.
- + **Easier to mill** geometry.
- + **Helical spiral** grooves for improved retention.
- + **Larger surface area for the prosthetic restoration** to avoid insufficient thickness depending on the materials used.
- + **Three different Ti-Base platforms:** Narrow, Regular, and Wide. The emergence profile will accordingly grow from implant to prosthetic restoration.
- + **Adjustable height:** With 8 mm as standard for the engaging versions and 6.5 mm for non-engaging.
- + **The recess for the Angled Screw Channel (ASC) can be adjusted** for both engaging and non-engaging. Check availability according to implant system.
- + **CAD-CAM** libraries are available for all possible variations.
- + Different **gingival heights** for up to 3.5 mm.
- + **Cutting guides** with fixing screw facilitate cutting and **trimming**



Be Precise


Scan Abutment

Accuracy, fit and reliability

One system, all solutions


Revolutionary high-precision Scan Abutment giving access to a fully integrated digital workflow.

A Scan Abutment suitable for both intraoral and desktop scanning devices enabling to perform implant level frameworks as well as giving access to our Custom Interface Ti-Base System or the one-of-a-kind high precision 3D printed models.



Scan Abutment

- + Compatible with intraoral systems (clinical) and with desktop systems (laboratory).
 - + Also compatible with the Renishaw® probe system.
 - + It includes a captive screw to prevent loss and make it easier to mount.
 - + The top section is made of PEEK, a polymer that favours implant reading and location.
 - + Maximum Z-precision, as it rests on the sagittal plane of the implant or analog.
 - + A tolerance of ± 5 microns throughout the entire manufacturing process.
 - + Reference engraved by laser etching.
 - + 2 heights for engaging: 10 mm for most cases, and 15 mm for situations involving a very deep implant.
- 1 Captive screw
 - 2 TiN coating
 - 3 Titanium connection
 - 4 Head PEEK
 - 5 Maximum Z-precision
 - 6 Reference laser etched
 - 7 2 heights for engaging versions



Optimal Best-fit

CAD Alignment

Digital accuracy means being able to align two STL meshes as accurately as possible.

IPD is the only manufacturer that has taken into account the oversizing that occurs when scanning and has developed a universal tool to help solve it.

By providing the customer with an assortment of 7 different STL files for each IPD scan bodies enables the possibility to improve CAD alignment regardless of the scanning device used.

Scan Transfer


The dedicated abutment for a precise IO-based workflow for full arch treatments. Scan Transfer is the first abutment offering clinicians a reliable protocol to precisely restore the most demanding clinical cases.

Its design and geometry allow them to be splinted creating a guide for the scanning accuracy by reducing scan distortion and increasing consistency and repeatability.

Thanks to Scan Transfer the most demanding clinical cases can be fully restored with an IO scanning device achieving the desired passive fit.



Technical Specifications



Verification Screw

Confirmation that the Scan Transfer is positioned correctly when the screw is flush with the top of the abutment.

TECHNICAL SPECIFICATIONS

- + Tapered design for easy axial scanning
- + 5 µm
- + Asymmetrical orientation cut-outs
- + Maximized scanning area
- + Retentive body that facilitates splinting

SPECIFIC LIBRARIES

- + Direct to implant library
- + Direct to implant library with ASC
- + Ti-Base level library
- + Ti-Base level library with ASC
- + Library for straight temporary abutment
- + Digital 3D Model

Compatible with:



Digital Analog


The most stable and precise

3D printing accuracy and reliability

The innovative design of the Digital Analog, with two fixation screws, ensures the exact positioning of the digital analog in a 3D printed model. The two fixation screws ensure that the analog is in the correct position, with one located laterally and the other positioned at the bottom of the 3D printed model. The lateral screw also allows the use of articulator without loss of the bite register.

Validation master

IPD provides a master STL file to allow final customer checking their printing device resolution. At the same time, a pretend 3D analog will be provided free of charge to allow testing the adjustment on the printed calibration master in a convenient, simple and fast way. Once acknowledge this value customer will have access to a customised CAD library with optimised parameter. There are existing up to 11 different CAD library files to ensure making it suitable to efficiently work with any 3D printing system.




Be Digital

CAD Solutions

For all the leading CAD software systems on the market.

Our aim is to simplify and improve the most your digital workflow by providing accurate solutions to your everyday digital challenges. We have developed and tested a full system to guarantee an accurate fully digital flow.

Scan Abutments Extra/Intra-Oral system is the definitive and most integrated solution on the market, allowing you to work with any type of CAD library, Ti-Bases, Direct to Implant (milled connection and sinter-milled connection), Digital Analog (3D printed model), Scan Transfer and more.



Ti-Base level full contour



Ti-Base level framework



Implant level sinter-milled connection



Implant level milled connection



Digital analog / 3D printed model

Implant Level Libraries

Library file that allows to freely design any framework type at implant level, with both straight and angled screw channel.



Ti-Base Level Libraries

Library file allowing the design of any type of prosthetic solution on top of the exclusive IPD Custom Ti-Base to produce cemented-screw retained restorations. It's a versatile solution for intended use with aesthetic raw materials offering the most appropriate final result for each clinical case thanks to its flexibility as well as all supported options. For example, the two adjustment tolerances: A tighter and a looser one to support from ceramic or PMMA frameworks for example.



3D digital model

Through intraoral or Desktop scanning will allow to create high precision 3D models simplifying and facilitating the whole clinical and lab protocol. The full system has been developed to provide a stable, accurate and cost-effective solution.





ModelPro is a dedicated CAD software to the design of 3D models, among many other options, offering in a single software the features of the main dental and mesh editing software's.

Cloud-based, it is the first CAD software that can run on both Windows and Mac OS®. The lack of a dongle allows unlimited installations, and the cloud-based environment offers free access via relevant credentials from anywhere, anytime. Exclusive for IPD customers it offers a top CAD software for designing 3D models completely free of charge.

Model Pro is also supporting unique features such as "Scan Abutment Replace", which allows switching to an IPD workflow from a third-party scan, or the exclusive "Stump Die 3D" system



- + 3D Model
- + Automatic detection margin prep from CAD
- + Mesh edit, reparation and restoration
- + Scan Abutment alignment
- + Detachable dies
- + Stump Die 3D
- + 3D articulator
- + Digital analog
- + Detachable soft tissue
- + Scan Abutment Replace



Index

Compatibilities

● AB Serie	MULTI-UNIT [ø4.8]	15
● NOBEL BIO CARE®		
AA Serie	BRANEMARK SYSTEM® [ø3.5 / ø4.1 / ø5.1]	18
AD Serie	NOBEL ACTIVE® [ø3.0 / ø3.5 / ø4.3]	21
● BIOMET® 3i®		
BA Serie	OSSEOTITE® [ø3.4 / ø4.1 / ø5]	24
BB Serie	CERTAIN® [ø3.4 / ø4.1 / ø5]	27
● STRAUMANN®		
DA Serie	TISSUE LEVEL [ø4.8 RN / ø6.5 WN]	30
DB Serie	BONE LEVEL [ø3.3 NC / ø4.1 RC]	33
DC Serie	BLX® [RB/WB]	36
● ASTRA®		
EA Serie	OSSEOSPEED™ [ø3 / ø3.5-4 / ø4.5-5]	39
EB Serie	EVOLUTION® [ø3.6 / ø4.2]	42
● ZIMMER®		
FA Serie	SCREW VENT® [ø3.5 / ø4.5 / ø5.7]	45
● BIOHORIZONS®		
LB Serie	TAPERED INTERNAL [ø3 / ø3.5 / ø4.5 / ø5.7]	48
● HIOSEN IMPLANT®		
OB Serie	TSIII [ø3.5 / ø4.0]	51
● NEODENT®		
RB Serie	GRAND MORSE™ GM [GM]	54
RC Serie	GM ABUTMENT [RP]	57
RD Serie	GM MICRO ABUTMENT [ø3.5]	59
● MIS®		
TA Serie	SEVEN® [Narrow 3.3 / Standard / Wide]	61
TB Serie	C1® [Standard]	64
● MEGAGEN®		
WA Serie	ANYONE® [RP]	67
WB Serie	ANYRIDGE® [ø4.0]	70

Analogs/Digital analogs

● NOBEL BIO CARE®		
AC Serie	NOBEL REPLACE® SELECT [ø3.5 / ø4.3 / ø5]	73
● KLOCKNER®		
CA Serie	SK2-NK2 [ø4.2]	73
CB Serie	ESSENTIAL CONE® [ø4.5]	73
CC Serie	KL™ [ø3.5 / ø4.1 / ø5.1]	74
CD Serie	VEGA® [ø3.5 / ø4.0 / ø4.5]	74
● STRAUMANN®		
DD Serie	SRA® [ø3.5/ø4.6]	74
● ZIMMER®		
FB Serie	SWISSPLUS® [ø4.8]	75
FC Serie	EZTETIC® [ø3.1]	75

●	BTI®		
	GA Serie	EXTERNAL [ø3.5 / ø4.1 / ø5.5]	75
	GB Serie	INTERNAL UNIVERSAL® [Universal® ø4.1 / Wide ø5.5]	76
	GC Serie	MULTI-IM® [ø4.1]	76
●	MICRODENT®		
	HA Serie	SYSTEM [ø3.5 / ø4.2 / ø5.1]	77
	HB Serie	UNIVERSAL™ [ø3.5 / ø4.1 / ø5.1]	77
●	DENTSPLY®		
	IA Serie	XIVE® FRIADENT® [ø3.4 / ø3.8 / ø4.5]	77
	IB Serie	ANKYLOS® [C/X]	77
●	CAMLOG®		
	JA Serie	CAMLOG® [ø3.3 / ø3.8 / ø4.3 / ø5.0]	78
	JB Serie	CONOLOG® [ø3.3 / ø3.8 / ø4.3]	78
●	SWEDEN & MARTINA®		
	MA Serie	OUTLINK® [ø3.3 / ø4.1 / ø5]	78
	MB Serie	PREMIUM™ KOHNO® [ø3.3 / ø3.8 / ø4.25 / ø5]	79
●	PHIBO®		
	NA Serie	TSH® [S2 / S3-S4 / S5]	79
●	BIOTECH® DENTAL		
	QB Serie	KONTACT® [RP]	80
●	NEODENT®		
	RA Serie	HELIX® HE [ø3.3 / ø4.1 / ø5.0]	80
●	BEGO® SEMADOS®		
	SB Serie	SC/RS [ø3.25 / ø3.75 / ø4.1 / ø4.5]	81
●	MIS®		
	TB Serie	V3® [Standard]	81
●	DIO®		
	UB Serie	UFII [NP / RP]	82
●	DENTIUM®		
	XA Serie	IMPLANTIUM® / SUPERLINE™ [RP]	82
●	MEDENTIS®		
	YB Serie	ICX [RP]	83
●	GLOBAL D®		
	1A Serie	IN-KONE® [RP]	83
●	ANTHOGYR®		
	2A Serie	AXIOM® BL [RP]	84

Working Systems

●	CUSTOM Ti-BASE SYSTEM	87
●	MULTI-UNIT SYSTEM	95
●	ACCESSORIES	101
●	TECHNICAL SPECIFICATIONS	103



Compatibilities



Compatibility

AB Serie MULTI-UNIT [ø4.8]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	ø4.8
	10	IPD/AB-SR-01

Captive screw included. 1.20 Hex Screwdriver.



SCAN TRANSFER

Indication	Height (mm)	Screw	ø4.8
	12.5	IPD/AB-TR-11/C	IPD/AB-SR-11

IPD/AB-SR-11 includes screw and a protective cap.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®

Indication	Post Height (mm)	Gingival Height	ø4.8
	3.5	0.2	IPD/AB-IR-01/3D
	4.0	0.2	IPD/AB-IR-11/3D
	6.5	0.2	IPD/AB-IR-03/3D

Please see page 87 for the Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	Height	ø4.8
-	Unigrip	Flat	-	IPD/AB-TR-00
			-	IPD/AB-TR-00/TIN
TiN	TPA	Flat	For 3.5 and 4.0 mm height ti-base	IPD/AB-TR-50
			For 6.5 mm height ti-base	IPD/AB-TR-51
	TPA	Flat	-	IPD/EA-TR-51

IPD/EA-TR-51 is the dedicated screw for MUA level restorations on Zr-PMMA.



ANALOGS

Types	ø4.8
Analog	IPD/AB-AR-00
Digital Analog	IPD/AB-AR-00/3D
Bottom Screw	IPD/KA-TA-00
Side Screw	IPD/KA-TA-01



OPEN AND CLOSED TRAY COPING

	ø4.8
Open tray coping	IPD/AB-CR-00
Closed tray coping	IPD/AB-CR-01



OVERDENTURE RESTORATIONS

PSD FOR MULTI-UNIT

Units	Height (mm)	ø4.8
2	3.35	IPD/AB-LR-02



PSD System Accessories

Indication	ø4.8
Analog	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00
Carrier	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT / OPEN TRAY COPING

Indication	ø4.8
	IPD/AB-PR-01
	IPD/AB-PR-00
	IPD/AB-PR-06
Coping screw	IPD/AB-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.
IPD/AB-PR-06 corresponds to the Classic design and cannot be used as a coping.



HEALING ABUTMENT

HEALING ABUTMENT

Material	ø4.8
Plastic POM	IPD/AB-DR-00
Titanium	IPD/AB-DR-01

IPD/AB-TR-00 screw is not included.



INSTRUMENTS AND ACCESSORIES

ANGLED MULTI-UNIT CARRIER

Indication	ø4.8
	IPD/LL-TM-00



SCREWDRIVERS

Type	Lengths (mm)	ø4.8
Unigrip	18	IPD/PU-00-18
	25	IPD/PU-00-25
	32	IPD/PU-00-32
2.0 Hex	Straight Multi-Unit abutment	20
		IPD/LL-AB-00



ACCESSORIES

Indication	Lengths (mm)	ø4.8
Clinical RA Shank Handle	16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90





Compatibility

NOBEL BIO CARE®

AA Serie BRANEMARK SYSTEM®

[$\phi 3.5$ / $\phi 4.1$ / $\phi 5.1$]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	$\phi 3.5$	$\phi 4.1$	$\phi 5.1$
	10	IPD/AA-SN-00	IPD/AA-SR-00	IPD/AA-SW-00
	15	IPD/AA-SN-02	IPD/AA-SR-02	-
	10	IPD/AA-SN-01	IPD/AA-SR-01	IPD/AA-SW-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®

Indication	Gingival Height	$\phi 3.5$	$\phi 4.1$	$\phi 5.1$
	0.5	-	IPD/AA-IR-00/3D	-
	1.5	IPD/AA-IN-02/3D	IPD/AA-IR-02/3D	-
	1.3	-	-	IPD/AA-IW-02/3D
	2.5	-	IPD/AA-IR-04/3D	-
	0.5	-	IPD/AA-IR-01/3D	-
	1.5	IPD/AA-IN-03/3D	IPD/AA-IR-03/3D	-
	1.3	-	-	IPD/AA-IW-03/3D
	2.5	-	IPD/AA-IR-05/3D	-


Please see page 87 for the Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	$\phi 3.5$	$\phi 4.1$	$\phi 5.1$
-	Unigrip	Conical	IPD/AA-TN-00	IPD/AA-TR-00	IPD/AA-TW-00
		Flat	-	IPD/AA-TR-02	-
	1.20 Hex	Low Profile conical	-	IPD/AA-TR-01	-
TiN	Unigrip	Conical	IPD/AA-TN-00/TIN	IPD/AA-TR-00/TIN	IPD/AA-TW-00/TIN
		Flat	-	IPD/AA-TR-02/TIN	-
	1.20 Hex	Low Profile conical	-	IPD/AA-TR-01/TIN	-
		Conical	IPD/AA-TN-50	IPD/AA-TR-50	IPD/AA-TW-50
	TPA	Flat	-	IPD/AA-TR-52	-

The IPD/AA-TW-01 screw is compatible with the Ti-Bases IPD/AA-IW-02/3D and IPD/AA-IW-03/3D



STRAIGHT MULTI-UNIT ABUTMENT

Angulation	Gingival Height	ø3.5	ø4.1	ø5.1
Straight	1.0	IPD/AA-MN-01	IPD/AB-MR-01	-
	2.0	IPD/AA-MN-02	IPD/AB-MR-02	-
	3.0	IPD/AA-MN-03	IPD/AB-MR-03	-
	4.0	-	IPD/AB-MR-04	-
	5.0	-	IPD/AB-MR-05	-

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	ø3.5	ø4.1	ø5.1
Analog	IPD/AA-AN-00	IPD/AA-AR-00	IPD/AA-AW-00
Digital Analog	IPD/AA-AN-00/3D	IPD/AA-AR-00/3D	IPD/AA-AW-00/3D
Bottom Screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Side screw	IPD/KA-TA-01	IPD/KA-TA-01	IPD/KA-TA-01



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	ø3.5	ø4.1	ø5.1
Analog	IPD/PS-AR-00	IPD/PS-AR-00	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Carrier	IPD/KA-CL-14	IPD/KA-CL-14	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	ø3.5	ø4.1	ø5.1
	IPD/AA-PN-06	IPD/AA-PR-06	IPD/AA-PW-06
	IPD/AA-PN-07	IPD/AA-PR-07	IPD/AA-PW-07
Coping screw	IPD/AA-TN-00/C	IPD/AA-TR-00/C	IPD/AA-TW-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	ø3.5	ø4.1	ø5.1
2.0	IPD/AA-DN-02	IPD/AA-DR-02	IPD/AA-DW-02
3.0	IPD/AA-DN-03	IPD/AA-DR-03	IPD/AA-DW-03
4.0	IPD/AA-DN-04	IPD/AA-DR-04	IPD/AA-DW-04
5.0	-	IPD/AA-DR-05	-



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	ø3.5	ø4.1	ø5.1
Unigrip	18	IPD/PU-00-18	IPD/PU-00-18	IPD/PU-00-18
	25	IPD/PU-00-25	IPD/PU-00-25	IPD/PU-00-25
	32	IPD/PU-00-32	IPD/PU-00-32	IPD/PU-00-32
1.20 Hex	18	IPD/PH-20-18	IPD/PH-20-18	IPD/PH-20-18
	25	IPD/PH-20-25	IPD/PH-20-25	IPD/PH-20-25
	32	IPD/PH-20-32	IPD/PH-20-32	IPD/PH-20-32
TPA	18	IPD/KA-CT-18	IPD/KA-CT-18	IPD/KA-CT-18
	25	IPD/KA-CT-25	IPD/KA-CT-25	IPD/KA-CT-25
	32	IPD/KA-CT-32	IPD/KA-CT-32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	ø3.5	ø4.1	ø5.1
Clinical RA Shank Handle	16	IPD/MC-00-16	IPD/MC-00-16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70	IPD/ML-00-70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90	IPD/DI-00-90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90	IPD/DI-KA-90	IPD/DI-KA-90





Compatibility

NOBEL BIO CARE®

AD Serie NOBEL ACTIVE®

[$\phi 3.0$ / $\phi 3.5$ / $\phi 4.3$]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	$\phi 3.0$	$\phi 3.5$	$\phi 4.3$
	10	IPD/AD-ST-00	IPD/AD-SN-00	IPD/AD-SR-00
	15	IPD/AD-ST-02	IPD/AD-SN-02	IPD/AD-SR-02
	10	IPD/AD-ST-01	IPD/AD-SN-01	IPD/AD-SR-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®

Indication	Gingival Height	$\phi 3.0$	$\phi 3.5$	$\phi 4.3$
	0.5	-	IPD/AD-IN-00/3D	IPD/AD-IR-00/3D
	0.8	IPD/AD-IT-00/3D	-	-
	1.5	IPD/AD-IT-02/3D	IPD/AD-IN-02/3D	IPD/AD-IR-02/3D
	2.5	IPD/AD-IT-04/3D	IPD/AD-IN-04/3D	IPD/AD-IR-04/3D
	3.5	-	IPD/AD-IN-06/3D	-
	0.5	-	IPD/AD-IN-01/3D	IPD/AD-IR-01/3D
	0.8	IPD/AD-IT-01/3D	-	-
	1.5	IPD/AD-IT-03/3D	IPD/AD-IN-03/3D	IPD/AD-IR-03/3D
	2.5	IPD/AD-IT-05/3D	IPD/AD-IN-05/3D	IPD/AD-IR-05/3D
	3.5	-	IPD/AD-IN-07/3D	-

High version of the non-engaging Custom Ti-Base subject to availability.

Please see page 87 for the Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	$\phi 3.0$	$\phi 3.5$	$\phi 4.3$
-	Unigrip	Conical	IPD/AD-TT-00	IPD/AA-TN-00	IPD/AA-TR-00
TiN	Unigrip	Conical	IPD/AD-TT-00/TIN	IPD/AA-TN-00/TIN	IPD/AA-TR-00/TIN
	TPA	Conical	IPD/AD-TT-50	IPD/AA-TN-50	IPD/AD-TR-50



STRAIGHT AND ANGLED MULTI-UNIT ABUTMENT

Angulation	Gingival Height	ø3.0	ø3.5	ø4.3
Straight	1.5	-	IPD/AD-MN-01	IPD/AD-MR-01
	2.5	-	IPD/AD-MN-02	IPD/AD-MR-02
	3.5	-	IPD/AD-MN-03	IPD/AD-MR-03
	4.5	-	IPD/AD-MN-04	IPD/AD-MR-04
17°	2.5	-	IPD/AD-MN-12	IPD/AD-MR-12
	3.5	-	IPD/AD-MN-13	IPD/AD-MR-13
30°	3.5	-	IPD/AD-MN-33	IPD/AD-MR-33
	4.5	-	IPD/AD-MN-34	IPD/AD-MR-34

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	ø3.0	ø3.5	ø4.3
Analog	IPD/AD-AT-00	IPD/AD-AN-00	IPD/AD-AR-00
Digital Analog	IPD/AD-AT-00/3D	IPD/AD-AN-00/3D	IPD/AD-AR-00/3D
Bottom Screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Side screw	IPD/KA-TA-01	IPD/KA-TA-01	IPD/KA-TA-01



OPEN TRAY COPING

Indication	ø3.0	ø3.5	ø4.3
Open tray coping	IPD/AD-CT-00	-	-

Classic design.



OVERDENTURE RESTORATIONS

PSD SYSTEM

Indication	Height (mm)	ø3.0	ø3.5	ø4.3
	1.0	-	IPD/AD-LN-01	IPD/AD-LR-01
	2.0	-	IPD/AD-LN-02	IPD/AD-LR-02
	3.0	-	IPD/AD-LN-03	IPD/AD-LR-03
	4.0	-	IPD/AD-LN-04	IPD/AD-LR-04
	5.0	-	IPD/AD-LN-05	IPD/AD-LR-05



PSD System Accessories

Indication	ø3.0	ø3.5	ø4.3
Analog	IPD/PS-AR-00	IPD/PS-AR-00	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Carrier	IPD/KA-CL-14	IPD/KA-CL-14	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	ø3.0	ø3.5	ø4.3
	IPD/AD-PT-06	IPD/AD-PN-06	IPD/AD-PR-06
	IPD/AD-PT-07	IPD/AD-PN-07	IPD/AD-PR-07
Coping screw	IPD/AD-TT-00/C	IPD/AD-TN-00/C	IPD/AD-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	ø3.0	ø3.5	ø4.3
1.5	-	IPD/AD-DN-01	IPD/AD-DR-01
3.0	IPD/AD-DT-03	IPD/AD-DN-03	IPD/AD-DR-03
5.0	IPD/AD-DT-05	IPD/AD-DN-05	IPD/AD-DR-05
7.0	IPD/AD-DT-07	-	-
3.0 wide	-	IPD/AD-DN-13	-
5.0 wide	-	IPD/AD-DN-15	-
Sealing cap	-	IPD/AD-DN-00	IPD/AD-DR-00



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	ø3.5	ø4.1	ø5.1
Unigrip	18	IPD/PU-00-18	IPD/PU-00-18	IPD/PU-00-18
	25	IPD/PU-00-25	IPD/PU-00-25	IPD/PU-00-25
	32	IPD/PU-00-32	IPD/PU-00-32	IPD/PU-00-32
1.20 Hex	18	IPD/PH-20-18	IPD/PH-20-18	IPD/PH-20-18
	25	IPD/PH-20-25	IPD/PH-20-25	IPD/PH-20-25
	32	IPD/PH-20-32	IPD/PH-20-32	IPD/PH-20-32
TPA	18	IPD/KA-CT-18	IPD/KA-CT-18	IPD/KA-CT-18
	25	IPD/KA-CT-25	IPD/KA-CT-25	IPD/KA-CT-25
	32	IPD/KA-CT-32	IPD/KA-CT-32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	ø3.5	ø4.1	ø5.1
Clinical RA Shank Handle	16	IPD/MC-00-16	IPD/MC-00-16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70	IPD/ML-00-70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90	IPD/DI-00-90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90	IPD/DI-KA-90	IPD/DI-KA-90





Compatibility

BIOMET® 3i®

BA Serie OSSEOTITE®

[$\phi 3.4$ / $\phi 4.1$ / $\phi 5.0$]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	$\phi 3.4$	$\phi 4.1$	$\phi 5.0$
	10	IPD/BA-SN-00	IPD/AA-SR-00	IPD/BA-SW-00
	15	IPD/BA-SN-02	IPD/AA-SR-02	IPD/BA-SW-02
	10	IPD/BA-SN-01	IPD/AA-SR-01	IPD/BA-SW-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®


Indication	Gingival Height	$\phi 3.4$	$\phi 4.1$	$\phi 5.0$
	0.5	-	IPD/AA-IR-00/3D	IPD/BA-IW-00/3D
	1.5	IPD/BA-IN-02/3D	IPD/AA-IR-02/3D	IPD/BA-IW-02/3D
	2.5	-	IPD/AA-IR-04/3D	-
	0.5	-	IPD/AA-IR-01/3D	IPD/BA-IW-01/3D
	1.5	IPD/BA-IN-03/3D	PD/AA-IR-03/3D	IPD/BA-IW-03/3D
	2.5	-	IPD/AA-IR-05/3D	-

High version of the non-engaging Custom Ti-Base subject to availability.
Please see page 87 for the Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	$\phi 3.4$	$\phi 4.1$	$\phi 5.0$
-	1.20 Hex	Conical	IPD/BA-TR-00	IPD/BA-TR-00	IPD/BA-TR-00
	1.25 Hex	Conical	IPD/BA-TR-01	IPD/BA-TR-01	IPD/BA-TR-01
	Square	Flat	IPD/BA-TR-02	IPD/BA-TR-02	IPD/BA-TR-02
	1.20 Hex	Low profile conical	IPD/AA-TR-01	IPD/AA-TR-01	IPD/AA-TR-01
TiN	1.20 Hex	Conical	IPD/BA-TR-00/TiN	IPD/BA-TR-00/TiN	IPD/BA-TR-00/TiN
	1.25 Hex	Conical	IPD/BA-TR-01/TiN	IPD/BA-TR-01/TiN	IPD/BA-TR-01/TiN
	Square	Flat	IPD/BA-TR-02/TiN	IPD/BA-TR-02/TiN	IPD/BA-TR-02/TiN
	1.20 Hex	Low profile conical	IPD/AA-TR-01/TiN	IPD/AA-TR-01/TiN	IPD/AA-TR-01/TiN
	TPA	Conical	IPD/AA-TR-50	IPD/AA-TR-50	IPD/AA-TR-50
		Flat	IPD/AA-TR-52	IPD/AA-TR-52	IPD/AA-TR-52



STRAIGHT MULTI-UNIT ABUTMENT

Angulation	Gingival Height	ø3.4	ø4.1	ø5.0
Straight	1.0	-	IPD/AB-MR-01	-
	2.0	-	IPD/AB-MR-02	-
	3.0	-	IPD/AB-MR-03	-
	4.0	-	IPD/AB-MR-04	-
	5.0	-	IPD/AB-MR-05	-

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	ø3.4	ø4.1	ø5.0
Analog	IPD/BA-AN-00	IPD/AA-AR-00	IPD/BA-AW-00
Digital Analog	IPD/BA-AN-00/3D	IPD/AA-AR-00/3D	IPD/BA-AW-00/3D
Bottom Screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Side screw	IPD/KA-TA-01	IPD/KA-TA-01	IPD/KA-TA-01



OPEN AND CLOSED TRAY COPING

	ø3.4	ø4.1	ø5.0
Open tray coping	-	IPD/AA-CR-00	IPD/BA-CW-00
Closed tray coping	-	IPD/AA-CR-01	-

Classic design.



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	ø3.4	ø4.1	ø5.0
Analog	IPD/PS-AR-00	IPD/PS-AR-00	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Carrier	IPD/KA-CL-14	IPD/KA-CL-14	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT / OPEN TRAY COPING

Indication	ø3.4	ø4.1	ø5.0
	IPD/BA-PN-06	IPD/AA-PR-06	-
	IPD/BA-PN-07	IPD/AA-PR-07	-
Coping screw	IPD/AA-TR-00/C	IPD/AA-TR-00/C	-

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	ø3.4	ø4.1	ø5.0
2.0	IPD/BA-DN-02	IPD/AA-DR-02	IPD/BA-DW-02
3.0	IPD/BA-DN-03	IPD/AA-DR-03	IPD/BA-DW-03
4.0	IPD/BA-DN-04	IPD/AA-DR-04	IPD/BA-DW-04
5.0	-	IPD/AA-DR-05	-



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	ø3.4	ø4.1	ø5.0
1.20 Hex	18	IPD/PH-20-18	IPD/PH-20-18	IPD/PH-20-18
	25	IPD/PH-20-25	IPD/PH-20-25	IPD/PH-20-25
	32	IPD/PH-20-32	IPD/PH-20-32	IPD/PH-20-32
1.25 Hex	18	IPD/PH-25-18	IPD/PH-25-18	IPD/PH-25-18
	25	IPD/PH-25-25	IPD/PH-25-25	IPD/PH-25-25
	32	IPD/PH-25-32	IPD/PH-25-32	IPD/PH-25-32
1.25 Square	18	IPD/PC-25-18	IPD/PC-25-18	IPD/PC-25-18
	25	IPD/PC-25-25	IPD/PC-25-25	IPD/PC-25-25
	32	IPD/PC-25-32	IPD/PC-25-32	IPD/PC-25-32
TPA	18	IPD/KA-CT-18	IPD/KA-CT-18	IPD/KA-CT-18
	25	IPD/KA-CT-25	IPD/KA-CT-25	IPD/KA-CT-25
	32	IPD/KA-CT-32	IPD/KA-CT-32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	ø3.4	ø4.1	ø5.0
Clinical RA Shank Handle	16	IPD/MC-00-16	IPD/MC-00-16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70	IPD/ML-00-70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90	IPD/DI-00-90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90	IPD/DI-KA-90	IPD/DI-KA-90





Compatibility

BIOMET® 3i®

BB Serie Certain®

[ø3.4 / ø4.1 / ø5.0]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	ø3.4	ø4.1	ø5.0
	10	IPD/BB-SN-00	IPD/BB-SR-00	IPD/BB-SW-00
	15	IPD/BB-SN-02	IPD/BB-SR-02	IPD/BB-SW-02
	10	IPD/BB-SN-01	IPD/BB-SR-01	IPD/BB-SW-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®


Indication	Gingival Height	ø3.4	ø4.1	ø5.0
	0.5	IPD/BB-IN-00/3D	IPD/BB-IR-00/3D	IPD/BB-IW-00/3D
	1.5	IPD/BB-IN-02/3D	IPD/BB-IR-02/3D	IPD/BB-IW-02/3D
	2.5	IPD/BB-IN-04/3D	IPD/BB-IR-04/3D	-
	0.5	IPD/BB-IN-01/3D	IPD/BB-IR-01/3D	IPD/BB-IW-01/3D
	1.5	IPD/BB-IN-03/3D	IPD/BB-IR-03/3D	IPD/BB-IW-03/3D
	2.5	IPD/BB-IN-05/3D	IPD/BB-IR-05/3D	-

High version of the non-engaging Custom Ti-Base subject to availability.
Please see page 87 for the  Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	ø3.4	ø4.1	ø5.0
-	1.20 Hex	Flat	IPD/BB-TR-01	IPD/BB-TR-01	IPD/BB-TR-01
		Conical	-	IPD/BB-TR-02	IPD/BB-TR-02
TiN	1.20 Hex	Flat	IPD/BB-TR-01/TIN	IPD/BB-TR-01/TIN	IPD/BB-TR-01/TIN
	TPA		IPD/BB-TR-51	IPD/BB-TR-51	IPD/BB-TR-51



STRAIGHT MULTI-UNIT ABUTMENT

Angulation	Height (mm)	ø3.4	ø4.1	ø5.0
Straight	1.0	IPD/BB-MN-01	IPD/BB-MR-01	-
	2.0	IPD/BB-MN-02	IPD/BB-MR-02	-
	3.0	IPD/BB-MN-03	IPD/BB-MR-03	-
	4.0	IPD/BB-MN-04	IPD/BB-MR-04	-
	5.0	-	IPD/BB-MR-05	-

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	ø3.4	ø4.1	ø5.0
Analog	IPD/BB-AN-00	IPD/BB-AR-00	IPD/BB-AW-00
Digital Analog	IPD/BB-AR-00/3D	IPD/BB-AR-00/3D	IPD/BB-AW-00/3D
Bottom Screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Side screw	IPD/KA-TA-01	IPD/KA-TA-01	IPD/KA-TA-01



OPEN TRAY COPING

	ø3.4	ø4.1	ø5.0
Open tray coping	-	-	IPD/BB-CW-00

Classic design.



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	ø3.4	ø4.1	ø5.0
Analog	IPD/PS-AR-00	IPD/PS-AR-00	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Carrier	IPD/KA-TA-01	IPD/KA-TA-01	IPD/KA-TA-01



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	ø3.4	ø4.1	ø5.0
	IPD/BB-PN-06	IPD/BB-PR-06	-
	IPD/BB-PN-07	IPD/BB-PR-07	-
Coping screw	IPD/BB-TR-01/C	IPD/BB-TR-01/C	-

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	ø3.4	ø4.1	ø5.0
2.0	IPD/BB-DN-02	IPD/BB-DR-02	IPD/BB-DW-02
3.0	IPD/BB-DN-03	IPD/BB-DR-03	IPD/BB-DW-03
4.0	IPD/BB-DN-04	IPD/BB-DR-04	IPD/BB-DW-04
5.0	-	IPD/BB-DR-05	-



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	ø3.4	ø4.1	ø5.0
1.20 Hex	18	IPD/PH-20-18	IPD/PH-20-18	IPD/PH-20-18
	25	IPD/PH-20-25	IPD/PH-20-25	IPD/PH-20-25
	32	IPD/PH-20-32	IPD/PH-20-32	IPD/PH-20-32
TPA	18	IPD/KA-CT-18	IPD/KA-CT-18	IPD/KA-CT-18
	25	IPD/KA-CT-25	IPD/KA-CT-25	IPD/KA-CT-25
	32	IPD/KA-CT-32	IPD/KA-CT-32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	ø3.4	ø4.1	ø5.0
Clinical RA Shank Handle	16	IPD/MC-00-16	IPD/MC-00-16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70	IPD/ML-00-70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90	IPD/DI-00-90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90	IPD/DI-KA-90	IPD/DI-KA-90





Compatibility

STRAUMANN®

DA Serie TISSUE LEVEL

[ø4.8 RN / ø6.5 WN]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	ø4.8 RN	ø6.5 WN
	10	IPD/DA-SR-00	IPD/DA-SW-00
	15	IPD/DA-SR-02	-
	10	IPD/DA-SR-01	IPD/DA-SW-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®

Indication	Gingival Height	ø4.8 RN	ø6.5 WN
	0.4	IPD/DA-IR-00/3D	IPD/DA-IW-02/3D
	0.4	IPD/DA-IR-01/3D	IPD/DA-IW-03/3D

Please see page 87 for the  Custom Ti Base System® accessories.




SCREW

Coating	Screwdriver	Seat	Compatibilities	ø4.8 RN	ø6.5 WN
-	Torx	Conical	Implant Level	IPD/DA-TR-01	IPD/DA-TR-01
		Flat	Implant Level	IPD/DA-TR-02	IPD/DA-TR-02
		Conical	SynOcta®	IPD/DA-TR-00	-
TiN	Torx	Conical	Implant Level	IPD/DA-TR-01/TIN	IPD/DA-TR-01/TIN
			Implant Level	IPD/DA-TR-03	IPD/DA-TR-03/TIN
			SynOcta®	IPD/DA-TR-00/TIN	-
	TPA		Implant Level	IPD/DA-TR-50	IPD/DA-TR-50
			SynOcta®	IPD/DA-TR-51	IPD/DA-TR-51

The IPD/DA-TR-03 screw has the same measurements as the original direct implant screw.

IPD/DA-TR-01 screw is a narrower headed screw to be used on  Custom Ti Base System®



STRAIGHT MULTI-UNIT ABUTMENT

Angulation	Gingival Height	ø4.8 RN	ø6.5 WN
Straight	0.9	IPD/DA-MR-01	-



ANALOGS

Types	Compatibilities	ø4.8 RN	ø6.5 WN
Analog		IPD/DA-AR-00	IPD/DA-AW-00
	SynOcta®	IPD/DA-AR-01	-
Digital Analog		IPD/DA-AR-00/3D	IPD/DA-AW-00/3D
	Bottom Screw	IPD/KA-TA-00	IPD/KA-TA-00
	Side screw	IPD/KA-TA-01	IPD/KA-TA-01



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	ø4.8 RN	ø6.5 WN
Analog	IPD/PS-AR-00	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D
	Bottom screw	IPD/KA-TA-00
	Carrier	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	Compatibilities	ø4.8 RN	ø6.5 WN
		IPD/DA-PR-06	IPD/DA-PW-06
	SynOcta®	IPD/DA-PR-09	-
		IPD/DA-PR-07	IPD/DA-PW-07
	SynOcta®	IPD/DA-PR-08	-
		IPD/DA-TR-01/C	IPD/DA-TR-01/C
Coping screw			

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	ø4.8 RN	ø6.5 WN
0.0 Sealing cap	IPD/DA-DR-00	-
1.5	IPD/DA-DR-01	-
3.0	IPD/DA-DR-03	-
4.5	IPD/DA-DR-04	-



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	ø4.8 RN	ø6.5 WN
1.20 Hex	18	IPD/PH-20-18	IPD/PH-20-18
	25	IPD/PH-20-25	IPD/PH-20-25
	32	IPD/PH-20-32	IPD/PH-20-32
TPA	18	IPD/KA-CT-18	IPD/KA-CT-18
	25	IPD/KA-CT-25	IPD/KA-CT-25
	32	IPD/KA-CT-32	IPD/KA-CT-32
Torx	18	IPD/PT-00-18	IPD/PT-00-18
	25	IPD/PT-00-25	IPD/PT-00-25
	32	IPD/PT-00-32	IPD/PT-00-32



ACCESSORIES

Indication	Lengths (mm)	ø4.8 RN	ø6.5 WN
Clinical RA Shank Handle	16	IPD/MC-00-16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90	IPD/DI-KA-90





Compatibility

STRAUMANN®

DB Serie BONE LEVEL

[ϕ 3.3 NC / ϕ 4.1 RC]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	ϕ 3.3 NC	ϕ 4.1 RC
	10	IPD/DB-SN-00	IPD/DB-SR-00
	15	IPD/DB-SN-02	IPD/DB-SR-02
	10	IPD/DB-SN-01	IPD/DB-SR-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE  Custom Interface System®

Indication	Height (mm)	ϕ 3.3 NC	ϕ 4.1 RC
	0.5	IPD/DB-IN-00/3D	IPD/DB-IR-00/3D
	1.5	IPD/DB-IN-02/3D	IPD/DB-IR-02/3D
	2.5	IPD/DB-IN-04/3D	IPD/DB-IR-04/3D
	3.5	IPD/DB-IN-06/3D	IPD/DB-IR-06/3D
	0.5	IPD/DB-IN-01/3D	IPD/DB-IR-01/3D
	1.5	IPD/DB-IN-03/3D	IPD/DB-IR-03/3D
	2.5	IPD/DB-IN-05/3D	IPD/DB-IR-05/3D
	3.5	IPD/DB-IN-07/3D	IPD/DB-IR-07/3D

Please see page 87 for the  Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	ϕ 3.3 NC	ϕ 4.1 RC
-	Torx	Conical	IPD/DB-TR-00	IPD/DB-TR-00
TiN	Torx	Conical	IPD/DB-TR-00/TiN	IPD/DB-TR-00/TiN
	TPA	Conical	IPD/DB-TR-50	IPD/DB-TR-50



STRAIGHT & ANGLED MULTI-UNIT ABUTMENT

Angulation	Gingival Height	ø3.3 NC	ø4.1 RC
Straight	1.5	IPD/DB-MN-01	IPD/DB-MR-01
	2.5	IPD/DB-MN-02	IPD/DB-MR-02
	3.5	IPD/DB-MN-03	IPD/DB-MR-03
	4.5	IPD/DB-MN-04	IPD/DB-MR-04
17°	2.5	IPD/DB-MN-12	IPD/DB-MR-12
	3.5	IPD/DB-MN-13	IPD/DB-MR-13
30°	3.5	IPD/DB-MN-33	IPD/DB-MR-33
	4.5	IPD/DB-MN-34	IPD/DB-MR-34

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	ø3.3 NC	ø4.1 RC
Analog	IPD/DB-AN-00	IPD/DB-AR-00
Digital Analog	IPD/DB-AN-00/3D	IPD/DB-AR-00/3D
Bottom Screw	IPD/KA-TA-00	IPD/KA-TA-00
Side screw	IPD/KA-TA-01	IPD/KA-TA-01



OVERDENTURE RESTORATIONS

PSD SYSTEM

Indication	Height (mm)	ø3.3 NC	ø4.1 RC
	1.0	IPD/DB-LN-01	IPD/DB-LR-01
	2.0	IPD/DB-LN-02	IPD/DB-LR-02
	3.0	IPD/DB-LN-03	IPD/DB-LR-03
	4.0	IPD/DB-LN-04	IPD/DB-LR-04
	5.0	IPD/DB-LN-05	IPD/DB-LR-05



PSD System Accessories

Indication	ø3.3 NC	ø4.1 RC
Analog	IPD/PS-AR-00	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00	IPD/KA-TA-00
Carrier	IPD/KA-CL-14	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	ø3.3 NC	ø4.1 RC
	IPD/DB-PN-06	IPD/DB-PR-06
	IPD/DB-PN-07	IPD/DB-PR-07
Coping screw	IPD/DB-TR-00/C	IPD/DB-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	ø3.3 NC	ø4.1 RC
2.0	IPD/DB-DN-02	IPD/DB-DR-02
3.5	IPD/DB-DN-03	-
4.0	-	IPD/DB-DR-04
5.0	IPD/DB-DN-05	-
6.0	-	IPD/DB-DR-06



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	ø3.3 NC	ø4.1 RC
Torx	18	IPD/PT-00-18	IPD/PT-00-18
	25	IPD/PT-00-25	IPD/PT-00-25
	32	IPD/PT-00-32	IPD/PT-00-32
1.20 Hex	18	IPD/PH-20-18	IPD/PH-20-18
	25	IPD/PH-20-25	IPD/PH-20-25
	32	IPD/PH-20-32	IPD/PH-20-32
TPA	18	IPD/KA-CT-18	IPD/KA-CT-18
	25	IPD/KA-CT-25	IPD/KA-CT-25
	32	IPD/KA-CT-32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	ø3.3 NC	ø4.1 RC
Clinical RA Shank Handle	16	IPD/MC-00-16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90	IPD/DI-KA-90





Compatibility

STRAUMANN®

DC Serie BLX®

[RB / WB]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	RB	WB
	10	IPD/DC-SR-00	IPD/DC-SR-00
	15	IPD/DC-SR-02	IPD/DC-SR-02
	10	IPD/DC-SR-01	IPD/DC-SR-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®

Indication	Gingival Height	RB	WB
	1.0	IPD/DC-IR-00/3D	IPD/DC-IR-00/3D
	2.0	IPD/DC-IR-02/3D	IPD/DC-IR-02/3D
	3.0	IPD/DC-IR-04/3D	IPD/DC-IR-04/3D
	1.0	IPD/DC-IR-01/3D	IPD/DC-IR-01/3D
	2.0	IPD/DC-IR-03/3D	IPD/DC-IR-03/3D
	3.0	IPD/DC-IR-05/3D	IPD/DC-IR-05/3D

Please see page 87 for the Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	RB	WB
-	Torx	Conical	IPD/DC-TR-00	IPD/DC-TR-00
TiN			IPD/DC-TR-01	IPD/DC-TR-01
			IPD/DC-TR-00/TiN	IPD/DC-TR-00/TiN

The IPD/DC-TR-01 model is an abutment extractor screw.



STRAIGHT & ANGLED MULTI-UNIT ABUTMENT

Angulation	Gingival Height	RB	WB
Straight	1.5	IPD/DC-MR-01	IPD/DC-MR-01
	2.5	IPD/DC-MR-02	IPD/DC-MR-02
	3.5	IPD/DC-MR-03	IPD/DC-MR-03
17°	3.5	IPD/DC-MR-13	IPD/DC-MR-13
30°	4.5	IPD/DC-MR-34	IPD/DC-MR-34

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	RB	WB
Analog	IPD/DC-AR-00	-
Digital Analog	IPD/DC-AR-00/3D	-
Bottom Screw	IPD/KA-TA-00	IPD/KA-TA-00
Side screw	IPD/KA-TA-01	IPD/KA-TA-01



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	RB	WB
Analog	IPD/PS-AR-00	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00	IPD/KA-TA-00
Carrier	IPD/KA-CL-14	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	RB	WB
	IPD/DC-PR-06	IPD/DC-PR-06
	IPD/DC-PR-07	IPD/DC-PR-07
Coping screw	IPD/DC-TR-00/C	IPD/DC-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	RB	WB
1.5	IPD/DC-DR-01	IPD/DC-DR-01
3.0	IPD/DC-DR-03	IPD/DC-DR-03
5.0	IPD/DC-DR-05	IPD/DC-DR-05



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	RB	WB
Torx	18	IPD/PT-00-18	IPD/PT-00-18
	25	IPD/PT-00-25	IPD/PT-00-25
	32	IPD/PT-00-32	IPD/PT-00-32
1.20 Hex	18	IPD/PH-20-18	IPD/PH-20-18
	25	IPD/PH-20-25	IPD/PH-20-25
	32	IPD/PH-20-32	IPD/PH-20-32



ACCESSORIES

Indication	Lengths (mm)	RB	WB
Clinical RA Shank Handle	16	IPD/MC-00-16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90	IPD/DI-KA-90





Compatibility

ASTRA[®]

EA Serie OSSEOSPEED[™]

[ϕ 3.0 / ϕ 3.5-4.0 / ϕ 4.5-5.0]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	ϕ 3.0	ϕ 3.5-4.0	ϕ 4.5-5.0
	10	IPD/EA-SN-00	IPD/EA-SR-00	IPD/EA-SW-00
	15	IPD/EA-SN-02	IPD/EA-SR-02	IPD/EA-SW-02
	10	IPD/EA-SN-01	IPD/EA-SR-01	IPD/EA-SW-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System[®]

Indication	Gingival Height	ϕ 3.0	ϕ 3.5-4.0	ϕ 4.5-5.0
	0.5	-	IPD/EA-IR-00/3D	IPD/EA-IW-00/3D
	0.8	IPD/EA-IN-00/3D	-	-
	1.5	IPD/EA-IN-02/3D	IPD/EA-IR-02/3D	IPD/EA-IW-02/3D
	2.5	IPD/EA-IN-04/3D	IPD/EA-IR-04/3D	IPD/EA-IW-04/3D
	3.5	-	IPD/EA-IR-06/3D	-
	0.5	-	IPD/EA-IR-01/3D	IPD/EA-IW-01/3D
	0.8	IPD/EA-IN-01/3D	-	-
	1.5	IPD/EA-IN-03/3D	IPD/EA-IR-03/3D	IPD/EA-IW-03/3D
	2.5	IPD/EA-IN-05/3D	IPD/EA-IR-05/3D	IPD/EA-IW-05/3D
	3.5	-	IPD/EA-IR-07/3D	-


High version of the non-engaging Custom Ti-Base subject to availability.
Please see page 87 for the  Custom Ti Base System[®] accessories.



SCREWS

Coating	Screwdriver	Seat	ø3.0	ø3.5-4.0	ø4.5-5.0
-	1.27 Hex	Conical	IPD/EA-TN-00	-	-
		Flat	-	IPD/EA-TR-00	IPD/EA-TW-00
		Flat	-	IPD/EA-TR-01	IPD/EA-TR-01
TiN	1.27 Hex	Conical	IPD/EA-TN-00/TiN	-	-
		Flat	-	IPD/EA-TR-00/TiN	IPD/EA-TW-00/TiN
		Flat	-	IPD/EA-TR-01/TiN	IPD/EA-TR-01/TiN
	TPA	Conical	IPD/EA-TN-50	-	-
		Flat	-	IPD/EA-TR-50	IPD/AC-TR-50
		Flat	-	-	IPD/AC-TR-51
		Flat	-	IPD/EA-TR-51	IPD/EA-TR-51

IPD/EA-TR-01 is a screw for the 20° and 45° UniAbutment.
IPD/AC-TR-50 is a specific screw for **IPD Custom Ti-Base System®** Ti-Bases.
IPD/EA-TR-51 is a screw for the 20° and 45° UniAbutment.



STRAIGHT MULTI-UNIT ABUTMENT

Angulation	Gingival Height	ø3.0	ø3.5-4.0	ø4.5-5.0
Straight	1.5	-	IPD/EA-MR-01	IPD/EA-MW-01
	2.5	-	IPD/EA-MR-02	IPD/EA-MW-02
	3.5	-	IPD/EA-MR-03	IPD/EA-MW-03
	4.5	-	IPD/EA-MR-04	IPD/EA-MW-04

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	ø3.0	ø3.5-4.0	ø4.5-5.0
Analog	IPD/EA-AN-00	IPD/EA-AR-00	IPD/EA-AW-00
	-	IPD/EA-AR-01	IPD/EA-AR-01
Digital Analog	IPD/EA-AN-00/3D	IPD/EA-AR-00/3D	IPD/EA-AW-00/3D
Bottom Screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Side screw	IPD/KA-TA-01	IPD/KA-TA-01	IPD/KA-TA-01

IPD/EA-AR-01 is an analog for 20° UniAbutment.



OPEN TRAY COPING

ø3.0	ø3.5-4.0	ø4.5-5.0
Open tray coping	IPD/EA-CN-00	-

Classic design.



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	ø3.0	ø3.5-4.0	ø4.5-5.0
Analog	IPD/PS-AR-00	IPD/PS-AR-00	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Carrier	IPD/KA-CL-14	IPD/KA-CL-14	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	ø3.0	ø3.5-4.0	ø4.5-5.0
	IPD/EA-PN-06	IPD/EA-PR-06	IPD/EA-PW-06
	IPD/EA-PN-07	IPD/EA-PR-07	IPD/EA-PW-07
Coping screw	IPD/EA-TN-00/C	IPD/EA-TR-00/C	IPD/EA-TW-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	ø3.0	ø3.5-4.0	ø4.5-5.0
2.0	-	IPD/EA-DR-02	IPD/EA-DW-02
4.0	-	IPD/EA-DR-04	IPD/EA-DW-04
6.0	-	IPD/EA-DR-06	IPD/EA-DW-06
8.0	-	IPD/EA-DR-08	IPD/EA-DW-08
10.0	-	IPD/EA-DR-10	-



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	ø3.0	ø3.5-4.0	ø4.5-5.0
1.27 Hex	18	IPD/PH-27-18	IPD/PH-27-18	IPD/PH-27-18
	25	IPD/PH-27-25	IPD/PH-27-25	IPD/PH-27-25
	32	IPD/PH-27-32	IPD/PH-27-32	IPD/PH-27-32
1.20 Hex	18	IPD/PH-20-18	IPD/PH-20-18	IPD/PH-20-18
	25	IPD/PH-20-25	IPD/PH-20-25	IPD/PH-20-25
	32	IPD/PH-20-32	IPD/PH-20-32	IPD/PH-20-32
TPA	18	IPD/KA-CT-18	IPD/KA-CT-18	IPD/KA-CT-18
	25	IPD/KA-CT-25	IPD/KA-CT-25	IPD/KA-CT-25
	32	IPD/KA-CT-32	IPD/KA-CT-32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	ø3.0	ø3.5-4.0	ø4.5-5.0
Clinical RA Shank Handle	16	IPD/MC-00-16	IPD/MC-00-16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70	IPD/ML-00-70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90	IPD/DI-00-90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90	IPD/DI-KA-90	IPD/DI-KA-90





Compatibility

ASTRA[®]

EB Serie EVOLUTION[®]

[$\phi 3.6$ / $\phi 4.2$]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	$\phi 3.6$	$\phi 4.2$
	10	IPD/EB-SN-00	IPD/EB-SR-00
	15	IPD/EB-SN-02	IPD/EB-SR-02
	10	IPD/EB-SN-01	IPD/EB-SR-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System[®]

Indication	Gingival Height	$\phi 3.6$	$\phi 4.2$
	0.5	IPD/EB-IN-00/3D	IPD/EB-IR-00/3D
	1.5	IPD/EB-IN-02/3D	IPD/EB-IR-02/3D
	2.5	IPD/EB-IN-04/3D	IPD/EB-IR-04/3D
	0.5	IPD/EB-IN-01/3D	IPD/EB-IR-01/3D
	1.5	IPD/EB-IN-03/3D	IPD/EB-IR-03/3D
	2.5	IPD/EB-IN-05/3D	IPD/EB-IR-05/3D

High version of the non-engaging Custom Ti-Base subject to availability.
Please see page 87 for the Custom Ti Base System[®] accessories.



SCREWS

Coating	Screwdriver	Seat	$\phi 3.6$	$\phi 4.2$
-	1.27 Hex	Conical	IPD/EB-TN-00	IPD/EB-TR-00
TiN	1.27 Hex	Conical	IPD/EB-TN-00/TiN	IPD/EB-TR-00/TiN
	TPA	Conical	IPD/EB-TN-50	IPD/EB-TR-50



STRAIGHT MULTI-UNIT ABUTMENT

Angulation	Gingival Height	ø3.6	ø4.2
Straight	1.5	IPD/EB-MN-01	IPD/EB-MR-01
	2.5	IPD/EB-MN-02	IPD/EB-MR-02
	3.5	IPD/EB-MN-03	IPD/EB-MR-03

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	ø3.6	ø4.2
Analog	IPD/EB-AN-00	IPD/EB-AR-00
Digital Analog	IPD/EB-AN-00/3D	IPD/EB-AR-00/3D
Bottom Screw	IPD/KA-TA-00	IPD/KA-TA-00
Side screw	IPD/KA-TA-01	IPD/KA-TA-01



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	ø3.6	ø4.2
Analog	IPD/PS-AR-00	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00	IPD/KA-TA-00
Carrier	IPD/KA-CL-14	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	ø3.6	ø4.2
	IPD/EB-PN-06	IPD/EB-PR-06
	IPD/EB-PN-07	IPD/EB-PR-07
Coping screw	IPD/EB-TN-00/C	IPD/EB-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	ø3.6	ø4.2
3.5	IPD/EB-DN-03	IPD/EB-DR-03
4.5	IPD/EB-DN-04	IPD/EB-DR-04
5.5	IPD/EB-DN-05	IPD/EB-DR-05



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	ø3.6	ø4.2
1.27 Hex	18	IPD/PH-27-18	IPD/PH-27-18
	25	IPD/PH-27-25	IPD/PH-27-25
	32	IPD/PH-27-32	IPD/PH-27-32
1.20 Hex	18	IPD/PH-20-18	IPD/PH-20-18
	25	IPD/PH-20-25	IPD/PH-20-25
	32	IPD/PH-20-32	IPD/PH-20-32
TPA	18	IPD/KA-CT-18	IPD/KA-CT-18
	25	IPD/KA-CT-25	IPD/KA-CT-25
	32	IPD/KA-CT-32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	ø3.6	ø4.2
Clinical RA Shank Handle	16	IPD/MC-00-16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90	IPD/DI-KA-90





Compatibility

ZIMMER®

FA Serie SCREW VENT®

[$\phi 3.5$ / $\phi 4.5$ / $\phi 5.7$]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	$\phi 3.5$	$\phi 4.5$	$\phi 5.7$
	10	IPD/FA-SN-00	IPD/FA-SR-00	IPD/FA-SW-00
	15	IPD/FA-SN-02	IPD/FA-SR-02	IPD/FA-SW-02
	10	IPD/FA-SN-01	IPD/FA-SR-01	IPD/FA-SW-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®

Indication	Gingival Height	$\phi 3.5$	$\phi 4.5$	$\phi 5.7$
	0.5	IPD/FA-IN-00/3D	IPD/FA-IR-00/3D	IPD/FA-IW-00/3D
	1.5	IPD/FA-IN-02/3D	IPD/FA-IR-02/3D	-
	2.5	IPD/FA-IN-04/3D	IPD/FA-IR-04/3D	-
	0.5	IPD/FA-IN-01/3D	IPD/FA-IR-01/3D	IPD/FA-IW-01/3D
	1.5	IPD/FA-IN-03/3D	IPD/FA-IR-03/3D	-
	2.5	IPD/FA-IN-05/3D	IPD/FA-IR-05/3D	-

Please see page 87 for the  Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	$\phi 3.5$	$\phi 4.5$	$\phi 5.7$
-	1.25 Hex	Conical	IPD/FA-TR-00	IPD/FA-TR-00	IPD/FA-TR-00
TiN	1.25 Hex	Conical	IPD/FA-TR-00/TiN	IPD/FA-TR-00/TiN	IPD/FA-TR-00/TiN
	TPA	Conical	IPD/FA-TR-50	IPD/FA-TR-50	IPD/FA-TR-50



STRAIGHT MULTI-UNIT ABUTMENT

Angulation	Gingival Height	ø3.5	ø4.5	ø5.7
Straight	0.9	IPD/FA-MN-00	-	-
	1.5	IPD/FA-MN-01	IPD/FA-MR-01	-
	2.5	IPD/FA-MN-02	IPD/FA-MR-02	-
	3.5	IPD/FA-MN-03	IPD/FA-MR-03	-
	4.5	IPD/FA-MN-04	IPD/FA-MR-04	-

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	ø3.5	ø4.5	ø5.7
Analog	IPD/FA-AN-00	IPD/FA-AR-00	IPD/FA-AW-00
Digital Analog	IPD/FA-AN-00/3D	IPD/FA-AR-00/3D	IPD/FA-AW-00/3D
Bottom Screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Side screw	IPD/KA-TA-01	IPD/KA-TA-01	IPD/KA-TA-01



OPEN TRAY COPING

	ø3.5	ø4.5	ø5.7
Open tray coping	IPD/FA-CN-00	-	IPD/FA-CW-00

Classic design.



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	ø3.5	ø4.5	ø5.7
Analog	IPD/PS-AR-00	IPD/PS-AR-00	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Carrier	IPD/KA-CL-14	IPD/KA-CL-14	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	ø3.5	ø4.5	ø5.7
	IPD/FA-PN-06	IPD/FA-PR-06	IPD/FA-PW-06
	IPD/FA-PN-07	IPD/FA-PR-07	IPD/FA-PW-07
Coping screw	IPD/FA-TR-00/C	IPD/FA-TR-00/C	IPD/FA-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	ø3.5	ø4.5	ø5.7
1.5	IPD/FA-DN-01	IPD/FA-DR-01	-
3.0	IPD/FA-DN-03	IPD/FA-DR-03	IPD/FA-DW-03
5.0	IPD/FA-DN-05	IPD/FA-DR-05	IPD/FA-DW-05



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	ø3.5	ø4.5	ø5.7
1.25 Hex	18	IPD/PH-25-18	IPD/PH-25-18	IPD/PH-25-18
	25	IPD/PH-25-25	IPD/PH-25-25	IPD/PH-25-25
	32	IPD/PH-25-32	IPD/PH-25-32	IPD/PH-25-32
1.20 Hex	18	IPD/PH-20-18	IPD/PH-20-18	IPD/PH-20-18
	25	IPD/PH-20-25	IPD/PH-20-25	IPD/PH-20-25
	32	IPD/PH-20-32	IPD/PH-20-32	IPD/PH-20-32
TPA	18	IPD/KA-CT-18	IPD/KA-CT-18	IPD/KA-CT-18
	25	IPD/KA-CT-25	IPD/KA-CT-25	IPD/KA-CT-25
	32	IPD/KA-CT-32	IPD/KA-CT-32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	ø3.5	ø4.5	ø5.7
Clinical RA Shank Handle	16	IPD/MC-00-16	IPD/MC-00-16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70	IPD/ML-00-70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90	IPD/DI-00-90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90	IPD/DI-KA-90	IPD/DI-KA-90





Compatibility

BIOHORIZONS®

LB Serie TAPERED INTERNAL

[ϕ 3.0 / ϕ 3.5 / ϕ 4.5 / ϕ 5.7]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	ϕ 3.0	ϕ 3.5	ϕ 4.5	ϕ 5.7
	10	IPD/LB-ST-00	IPD/LB-SN-00	IPD/LB-SR-00	IPD/LB-SW-00
	15	IPD/LB-ST-02	IPD/LB-SN-02	IPD/LB-SR-02	IPD/LB-SW-02
	10	IPD/LB-ST-01	IPD/LB-SN-01	IPD/LB-SR-01	IPD/LB-SW-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®

Indication	Gingival Height	ϕ 3.0	ϕ 3.5	ϕ 4.5	ϕ 5.7
	0.5	IPD/LB-IT-00/3D	IPD/LB-IN-00/3D	IPD/LB-IR-00/3D	IPD/LB-IW-00/3D
	1.5	IPD/LB-IT-02/3D	IPD/LB-IN-02/3D	IPD/LB-IR-02/3D	-
	2.5	IPD/LB-IT-04/3D	IPD/LB-IN-04/3D	IPD/LB-IR-04/3D	-
	0.5	IPD/LB-IT-01/3D	IPD/LB-IN-01/3D	IPD/LB-IR-01/3D	IPD/LB-IW-01/3D
	1.5	IPD/LB-IT-03/3D	IPD/LB-IN-03/3D	IPD/LB-IR-03/3D	-
	2.5	IPD/LB-IT-05/3D	IPD/LB-IN-05/3D	IPD/LB-IR-05/3D	-

Please see page 87 for the Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	ϕ 3.0	ϕ 3.5	ϕ 4.5	ϕ 5.7
-	1.27 Hex	Conical	IPD/LB-TR-00	IPD/LB-TR-00	IPD/LB-TR-00	IPD/LB-TR-00
TiN	1.27 Hex	Conical	IPD/LB-TR-00/TIN	IPD/LB-TR-00/TIN	IPD/LB-TR-00/TIN	IPD/LB-TR-00/TIN
	TPA	Conical	IPD/LB-TR-50	IPD/LB-TR-50	IPD/LB-TR-50	IPD/LB-TR-50



STRAIGHT & ANGLED MULTI-UNIT ABUTMENT

Angulation	Height (mm)	Ø3.0	Ø3.5	Ø4.5	Ø5.7
Straight	1.5	-	IPD/FA-MN-01	IPD/FA-MR-01	-
	2.5	-	IPD/FA-MN-02	IPD/FA-MR-02	-
	3.5	-	IPD/FA-MN-03	IPD/FA-MR-03	-
	4.5	-	IPD/FA-MN-04	IPD/FA-MR-04	-
17°	2.5	-	IPD/LB-MN-12	-	-
	3.5	-	IPD/LB-MN-13	-	-
30°	3.5	-	IPD/LB-MN-33	-	-
	4.5	-	IPD/LB-MN-34	-	-

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	Ø3.0	Ø3.5	Ø4.5	Ø5.7
Analog	IPD/LB-AT-00	IPD/LB-AN-00	IPD/LB-AR-00	IPD/LB-AW-00
Digital Analog	IPD/LB-AT-00/3D	IPD/LB-AN-00/3D	IPD/LB-AR-00/3D	IPD/LB-AW-00/3D
Bottom Screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Side screw	IPD/KA-TA-01	IPD/KA-TA-01	IPD/KA-TA-01	IPD/KA-TA-01



OPEN TRAY COPING

Indication	Ø3.0	Ø3.5	Ø4.5	Ø5.7
Open tray coping	-	-	-	IPD/LB-CW-00

Classic design.



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	Ø3.0	Ø3.5	Ø4.5	Ø5.7
Analog	IPD/PS-AR-00	IPD/PS-AR-00	IPD/PS-AR-00	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Carrier	IPD/KA-CL-14	IPD/KA-CL-14	IPD/KA-CL-14	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	Ø3.0	Ø3.5	Ø4.5	Ø5.7
	IPD/LB-PT-06	IPD/LB-PN-06	IPD/LB-PR-06	IPD/LB-PW-06
	IPD/LB-PT-07	IPD/LB-PN-07	IPD/LB-PR-07	IPD/LB-PW-07
Coping screw	IPD/LB-TR-00/C	IPD/LB-TR-00/C	IPD/LB-TR-00/C	IPD/LB-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	Ø3.0	Ø3.5	Ø4.5	Ø5.7
1.5	-	IPD/FA-DN-01	IPD/FA-DR-01	-
2.0	IPD/LB-DT-02	-	-	-
3.0	IPD/LB-DT-03	IPD/FA-DN-03	IPD/FA-DR-03	IPD/FA-DW-03
5.0	IPD/LB-DT-05	IPD/FA-DN-05	IPD/FA-DR-05	IPD/FA-DW-05



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	Ø3.0	Ø3.5	Ø4.5	Ø5.7
1.27 Hex	18	IPD/PH-27-18	IPD/PH-27-18	IPD/PH-27-18	IPD/PH-27-18
	25	IPD/PH-27-25	IPD/PH-27-25	IPD/PH-27-25	IPD/PH-27-25
	32	IPD/PH-27-32	IPD/PH-27-32	IPD/PH-27-32	IPD/PH-27-32
1.20 Hex	18	IPD/PH-20-18	IPD/PH-20-18	IPD/PH-20-18	IPD/PH-20-18
	25	IPD/PH-20-25	IPD/PH-20-25	IPD/PH-20-25	IPD/PH-20-25
	32	IPD/PH-20-32	IPD/PH-20-32	IPD/PH-20-32	IPD/PH-20-32
TPA	18	IPD/KA-CT-18	IPD/KA-CT-18	IPD/KA-CT-18	IPD/KA-CT-18
	25	IPD/KA-CT-25	IPD/KA-CT-25	IPD/KA-CT-25	IPD/KA-CT-25
	32	IPD/KA-CT-32	IPD/KA-CT-32	IPD/KA-CT-32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	Ø3.0	Ø3.5	Ø4.5	Ø5.7
Clinical RA Shank Handle	16	IPD/MC-00-16	IPD/MC-00-16	IPD/MC-00-16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70	IPD/ML-00-70	IPD/ML-00-70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90	IPD/DI-00-90	IPD/DI-00-90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90	IPD/DI-KA-90	IPD/DI-KA-90	IPD/DI-KA-90





Compatibility

HIOSSEN IMPLANT® (OSSTEM)

OB Serie TSIII

[Ø3.5 M / Ø4.0 R]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	Ø3.5 M	Ø4.0 R
	10	IPD/OB-SN-00	IPD/OB-SR-00
	15	IPD/OB-SN-02	IPD/OB-SR-02
	10	IPD/OB-SN-01	IPD/OB-SR-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®

Indication	Gingival Height	Ø3.5 M	Ø4.0 R
	1.0	IPD/OB-IN-00/3D	IPD/OB-IR-00/3D
	2.0	IPD/OB-IN-02/3D	IPD/OB-IR-02/3D
	3.0	IPD/OB-IN-04/3D	IPD/OB-IR-04/3D
	1.0	IPD/OB-IN-01/3D	IPD/OB-IR-01/3D
	2.0	IPD/OB-IN-03/3D	IPD/OB-IR-03/3D
	3.0	IPD/OB-IN-05/3D	IPD/OB-IR-05/3D

Please see page 87 for the Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	Ø3.5 M	Ø4.0 R
-	1.20 Hex	Conical	IPD/OB-TN-00	IPD/OB-TR-00
TiN	1.20 Hex	Conical	IPD/OB-TN-00/TiN	IPD/OB-TR-00/TiN
	TPA	Conical	IPD/OB-TN-50	IPD/OB-TR-50



STRAIGHT & ANGLED MULTI-UNIT ABUTMENT

Angulation	Gingival Height	Ø3.5 M	Ø4.0 R
Straight	1.5	IPD/OB-MN-01	IPD/OB-MR-01
	2.5	IPD/OB-MN-02	IPD/OB-MR-02
	3.5	IPD/OB-MN-03	IPD/OB-MR-03
	4.5	IPD/OB-MN-04	IPD/OB-MR-04
	5.5	IPD/OB-MN-05	IPD/OB-MR-05
17°	2.5	IPD/OB-MN-12	IPD/OB-MR-12
	3.5	IPD/OB-MN-13	IPD/OB-MR-13
30°	3.5	IPD/OB-MN-33	IPD/OB-MR-33
	4.5	IPD/OB-MN-34	IPD/OB-MR-34

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	Ø3.5 M	Ø4.0 R
Analog	IPD/OB-AN-00	IPD/OB-AR-00
Digital Analog	IPD/OB-AN-00/3D	IPD/OB-AR-00/3D
Bottom Screw	IPD/KA-TA-00	IPD/KA-TA-00
Side screw	IPD/KA-TA-01	IPD/KA-TA-01



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	Ø3.5 M	Ø4.0 R
Analog	IPD/PS-AR-00	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00	IPD/KA-TA-00
Carrier	IPD/KA-CL-14	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	Ø3.5 M	Ø4.0 R
	-	IPD/OB-PN-06
	-	IPD/OB-PN-07
Coping screw	IPD/OB-TN-00/C	IPD/OB-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	Ø3.5 M	Ø4.0 R
3.0	IPD/OB-DN-03	IPD/OB-DR-03
5.0	IPD/OB-DN-05	IPD/OB-DR-05
7.0	IPD/OB-DN-07	IPD/OB-DR-07



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	Ø3.5 M	Ø4.0 R
1.20 Hex	18	IPD/PH-20-18	IPD/PH-20-18
	25	IPD/PH-20-25	IPD/PH-20-25
	32	IPD/PH-20-32	IPD/PH-20-32
TPA	18	IPD/KA-CT-18	IPD/KA-CT-18
	25	IPD/KA-CT-25	IPD/KA-CT-25
	32	IPD/KA-CT-32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	Ø3.5 M	Ø4.0 R
Clinical RA Shank Handle	16	IPD/MC-00-16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90	IPD/DI-KA-90





Compatibility

NEODENT®

RB Serie GRAND MORSE™ GM
[GM]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	GM
	10	IPD/RB-SR-00
	15	IPD/RB-SR-02
	10	IPD/RB-SR-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®

Indication	Gingival Height	GM
	1.0	IPD/RB-IR-00/3D
	2.0	IPD/RB-IR-02/3D
	3.0	IPD/RB-IR-04/3D
	4.0	IPD/RB-IR-06/3D
	1.0	IPD/RB-IR-01/3D
	2.0	IPD/RB-IR-03/3D
	3.0	IPD/RB-IR-05/3D
	4.0	IPD/RB-IR-07/3D

Please see page 87 for the  Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	GM
-	Neo GM	Conical	IPD/RB-TR-00
TiN	Neo GM	Conical	IPD/RB-TR-00/TiN
	TPA	Conical	IPD/RB-TR-50



STRAIGHT & ANGLED MULTI-UNIT ABUTMENT

Angulation	Gingival Height	GM
Straight	1.5	IPD/RB-MR-01
	2.5	IPD/RB-MR-02
	3.5	IPD/RB-MR-03
	4.5	IPD/RB-MR-04
17°	2.5	IPD/RB-MR-12
30°	3.5	IPD/RB-MR-33

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	GM
Analog	IPD/RB-AR-00
Digital Analog	IPD/RB-AR-00/3D
Bottom Screw	IPD/KA-TA-00
Side screw	IPD/KA-TA-01



OVERDENTURE RESTORATIONS

PSD SYSTEM

Indication	Height (mm)	GM
	2.0	IPD/RB-LR-02
	3.0	IPD/RB-LR-03
	4.0	IPD/RB-LR-04
	5.0	IPD/RB-LR-05



PSD System Accessories

Indication	GM
Analog	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00
Carrier	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	GM
	IPD/RB-PR-06
	IPD/RB-PR-07
Coping screw	IPD/RB-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	GM
1.5	IPD/RB-DR-01
3.0	IPD/RB-DR-03
5.0	IPD/RB-DR-05
7.0	IPD/RB-DR-07



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	GM
1.20 Hex	18	IPD/PH-20-18
	25	IPD/PH-20-25
	32	IPD/PH-20-32
Neo GM	18	IPD/PU-01-18
	25	IPD/PU-01-25
	32	IPD/PU-01-32
TPA	18	IPD/KA-CT-18
	25	IPD/KA-CT-25
	32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	GM
Clinical RA Shank Handle	16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90
PSD and insert extractor tool	-	IPD/LL-PS-01





Compatibility

NEODENT®

RC Serie GM ABUTMENT

[RP]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	RP
	10	IPD/RC-SR-00

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®

Indication	Gingival Height	RP
	0.5	IPD/RC-IR-00/3D

Please see page 87 for the  Custom Ti Base System® accessories.
Use wide platform cutting guides.



SCREWS

Coating	Screwdriver	Seat	RP
-	Neo GM	Flat	IPD/RC-TR-00
TiN	Neo GM	Flat	IPD/RC-TR-00/TiN
	TPA	Flat	IPD/RC-TR-50



ANALOGS

Types	RP
Analog	IPD/RC-AR-00
Digital Analog	IPD/RC-AR-00/3D
Bottom Screw	IPD/KA-TA-00
Side screw	IPD/KA-TA-01



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	RP
Analog	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00
Carrier	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	RP
	IPD/RC-PR-06
Coping screw	IPD/RC-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	RP
4.75	IPD/RC-DR-04



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	RP
1.20 Hex	18	IPD/PH-20-18
	25	IPD/PH-20-25
	32	IPD/PH-20-32
Neo GM	18	IPD/PU-01-18
	25	IPD/PU-01-25
	32	IPD/PU-01-32
TPA	18	IPD/KA-CT-18
	25	IPD/KA-CT-25
	32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	RP
Clinical RA Shank Handle	16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90





Compatibility

NEODENT®

RD Serie GM MICRO ABUTMENT
[ø3.5]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	ø3.5
	10	IPD/RD-SN-00
	10	IPD/RD-SN-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®

Indication	Gingival Height	ø3.5
	0.5	IPD/RD-IN-00/3D
	0.3	IPD/RD-IN-01/3D

Please see page 87 for the  Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	ø3.5
-	Neo GM	Conical	IPD/RD-TN-00
TiN	Neo GM	Conical	IPD/RD-TN-00/TiN
	TPA	Conical	IPD/RD-TN-50



ANALOGS

Types	ø3.5
Analog	IPD/RD-AN-00
Digital Analog	IPD/RD-AN-00/3D
Bottom Screw	IPD/KA-TA-00
Side screw	IPD/KA-TA-01



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	ø3.5
Analog	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00
Carrier	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	ø3.5
	IPD/RD-PN-06
	IPD/RD-PN-07
Coping screw	IPD/RD-TN-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	ø3.5
3.25	IPD/RD-DN-03



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	ø3.5
1.20 Hex	18	IPD/PH-20-18
	25	IPD/PH-20-25
	32	IPD/PH-20-32
Neo GM	18	IPD/PU-01-18
	25	IPD/PU-01-25
	32	IPD/PU-01-32
TPA	18	IPD/KA-CT-18
	25	IPD/KA-CT-25
	32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	ø3.5
Clinical RA Shank Handle	16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90





Compatibility

MIS[®]

TA Serie SEVEN[®]

[Narrow 3.3 / Standard / Wide]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	Narrow 3.3	Standard	Wide
	10	IPD/TA-ST-00	IPD/FA-SN-00	IPD/FA-SR-00
	15	IPD/TA-ST-02	IPD/FA-SN-02	IPD/FA-SR-02
	10	IPD/TA-ST-01	IPD/FA-SN-01	IPD/FA-SR-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System[®]


Indication	Gingival Height	Narrow 3.3	Standard	Wide
	0.5	IPD/TA-IT-00/3D	IPD/FA-IN-00/3D	IPD/FA-IR-00/3D
	1.5	IPD/TA-IT-02/3D	IPD/FA-IN-02/3D	IPD/FA-IR-02/3D
	2.5	IPD/TA-IT-04/3D	IPD/FA-IN-04/3D	IPD/FA-IR-04/3D
	0.5	IPD/TA-IT-01/3D	IPD/FA-IN-01/3D	IPD/FA-IR-01/3D
	1.5	IPD/TA-IT-03/3D	IPD/FA-IN-03/3D	IPD/FA-IR-03/3D
	2.5	IPD/TA-IT-05/3D	IPD/FA-IN-05/3D	IPD/FA-IR-05/3D

Please see page 87 for the  Custom Ti Base System[®] accessories.



SCREWS

Coating	Screwdriver	Seat	Narrow 3.3	Standard	Wide
-	1.25 Hex	Conical	IPD/TA-TT-00	IPD/FA-TR-00	IPD/FA-TR-00
TiN	1.25 Hex	Conical	IPD/TA-TT-00/TiN	IPD/FA-TR-00/TiN	IPD/FA-TR-00/TiN
	TPA	Conical	IPD/TA-TT-50	IPD/FA-TR-50	IPD/FA-TR-50



STRAIGHT MULTI-UNIT ABUTMENT

Angulation	Gingival Height	Narrow 3.3	Standard	Wide
Straight	0.9	-	IPD/FA-MN-00	-
	1.5	-	IPD/FA-MN-01	IPD/FA-MR-01
	2.5	-	IPD/FA-MN-02	IPD/FA-MR-02
	3.5	-	IPD/FA-MN-03	IPD/FA-MR-03
	4.5	-	IPD/FA-MN-04	IPD/FA-MR-04

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	Narrow 3.3	Standard	Wide
Analog	IPD/TA-AT-00	IPD/FA-AN-00	IPD/FA-AR-00
Digital Analog	IPD/TA-AT-00/3D	IPD/FA-AN-00/3D	IPD/FA-AR-00/3D
Bottom Screw	IPD/KA-TA-00	-	-
Side screw	IPD/KA-TA-01	-	-



OPEN TRAY COPING

	Narrow 3.3	Standard	Wide
Open tray coping	-	IPD/FA-CN-00	-

Classic design.



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	Narrow 3.3	Standard	Wide
Analog	IPD/PS-AR-00	IPD/PS-AR-00	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00	IPD/KA-TA-00	IPD/KA-TA-00
Carrier	IPD/KA-CL-14	IPD/KA-CL-14	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	Narrow 3.3	Standard	Wide
	IPD/TA-PT-06	IPD/FA-PN-06	IPD/FA-PR-06
	IPD/TA-PT-07	IPD/FA-PN-07	IPD/FA-PR-07
Coping screw	IPD/TA-TT-00/C	IPD/FA-TR-00/C	IPD/FA-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	Narrow 3.3	Standard	Wide
1.5	-	IPD/FA-DN-01	IPD/FA-DR-01
2.0	IPD/TA-DT-02	-	-
3.0	IPD/TA-DT-03	IPD/FA-DN-03	IPD/FA-DR-03
4.0	IPD/TA-DT-04	-	-
5.0	IPD/TA-DT-05	IPD/FA-DN-05	IPD/FA-DR-05



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	Narrow 3.3	Standard	Wide
1.25 Hex	18	IPD/PH-25-18	IPD/PH-25-18	IPD/PH-25-18
	25	IPD/PH-25-25	IPD/PH-25-25	IPD/PH-25-25
	32	IPD/PH-25-32	IPD/PH-25-32	IPD/PH-25-32
1.20 Hex	18	IPD/PH-20-18	IPD/PH-20-18	IPD/PH-20-18
	25	IPD/PH-20-25	IPD/PH-20-25	IPD/PH-20-25
	32	IPD/PH-20-32	IPD/PH-20-32	IPD/PH-20-32
TPA	18	IPD/KA-CT-18	IPD/KA-CT-18	IPD/KA-CT-18
	25	IPD/KA-CT-25	IPD/KA-CT-25	IPD/KA-CT-25
	32	IPD/KA-CT-32	IPD/KA-CT-32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	Narrow 3.3	Standard	Wide
Clinical RA Shank Handle	16	IPD/MC-00-16	IPD/MC-00-16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70	IPD/ML-00-70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90	IPD/DI-00-90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90	IPD/DI-KA-90	IPD/DI-KA-90





Compatibility

MIS[®]

TB Serie C1[®]

[Standard]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	Standard
	10	IPD/TB-SR-00
	15	IPD/TB-SR-02
	10	IPD/TB-SR-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System[®]

Indication	Gingival Height	Standard
	0.7	IPD/TB-IR-00/3D
	1.5	IPD/TB-IR-02/3D
	2.5	IPD/TB-IR-04/3D
	3.5	IPD/TB-IR-06/3D
	0.7	IPD/TB-IR-01/3D
	1.5	IPD/TB-IR-03/3D
	2.5	IPD/TB-IR-05/3D
	3.5	IPD/TB-IR-07/3D

Please see page 87 for the  Custom Ti Base System[®] accessories.
High version of the non-engaging Custom Ti-Base subject to availability.



SCREWS

Coating	Screwdriver	Seat	Standard
-	1.25 Hex	Conical	IPD/TB-TR-00
	1.25 Hex	Conical	IPD/TB-TR-00/TIN
TiN	TPA	Conical	IPD/TB-TR-50



STRAIGHT MULTI-UNIT ABUTMENT

Angulation	Gingival Height	Standard
Straight	1.0	IPD/TB-MR-01
	2.0	IPD/TB-MR-02
	3.0	IPD/TB-MR-03
	4.0	IPD/TB-MR-04

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	Standard
Analog	IPD/TB-AR-00
Digital Analog	IPD/TB-AR-00/3D
Bottom Screw	IPD/KA-TA-00
Side screw	IPD/KA-TA-01



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	Standard
Analog	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00
Carrier	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	Standard
	IPD/TB-PR-06
	IPD/TB-PR-07
Coping screw	IPD/TB-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	Standard
2.0	IPD/TB-DR-02
3.0	IPD/TB-DR-03
4.0	IPD/TB-DR-04
5.0	IPD/TB-DR-05



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	Standard
1.25 Hex	18	IPD / PH-25-18
	25	IPD / PH-25-25
	32	IPD / PH-25-32
1.20 Hex	18	IPD / PH-20-18
	25	IPD / PH-20-25
	32	IPD / PH-20-32
TPA	18	IPD / KA-CT-18
	25	IPD / KA-CT-25
	32	IPD / KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	Standard
Clinical RA Shank Handle	16	IPD / MC-00-16
Laboratory RA Shank Handle	70	IPD / ML-00-70
RA Shank Ratchet	90	IPD / DI-00-90
Ratchet adapter	-	IPD / DI-KA-90





Compatibility

MEGAGEN®

WA Serie ANYONE®

[RP]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	RP
	10	IPD/WA-SR-00
	15	IPD/WA-SR-02
	10	IPD/WA-SR-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®

Indication	Gingival Height	RP
	1.0	IPD/WA-IR-00/3D
	2.0	IPD/WA-IR-02/3D
	3.0	IPD/WA-IR-04/3D
	1.0	IPD/WA-IR-01/3D
	2.0	IPD/WA-IR-03/3D
	3.0	IPD/WA-IR-05/3D

Please see page 87 for the  Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	RP
-	1.20 Hex	Conical	IPD/WA-TR-00
TiN	1.20 Hex	Conical	IPD/WA-TR-00/TiN
	TPA	Conical	IPD/WA-TR-50



STRAIGHT MULTI-UNIT ABUTMENT

Angulation	Height (mm)	RP
Straight	1.5	IPD/WA-MR-01
	2.5	IPD/WA-MR-02
	3.5	IPD/WA-MR-03

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	RP
Analog	IPD/WA-AR-00
Digital Analog	IPD/WA-AR-00/3D
Bottom Screw	IPD/KA-TA-00
Side screw	IPD/KA-TA-01



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	RP
Analog	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00
Carrier	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	RP
	IPD/WA-PR-06
	IPD/WA-PR-07
Coping screw	IPD/WA-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	RP
3.0	IPD/WA-DR-03
5.0	IPD/WA-DR-05
7.0	IPD/WA-DR-07



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS

Type	Lengths (mm)	RP
1.20 Hex	18	IPD/PH-20-18
	25	IPD/PH-20-25
	32	IPD/PH-20-32
TPA	18	IPD/KA-CT-18
	25	IPD/KA-CT-25
	32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	RP
Clinical RA Shank Handle	16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90





Compatibility

MEGAGEN®

WB Serie ANYRIDGE®

[ø4.0]



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	ø4.0
	10	IPD/WB-SR-00
	15	IPD/WB-SR-02
	10	IPD/WB-SR-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE Custom Ti-Base System®

Indication	Gingival Height	ø4.0
	1.0	IPD/WB-IR-00/3D
	2.0	IPD/WB-IR-02/3D
	3.0	IPD/WB-IR-04/3D
	1.0	IPD/WB-IR-01/3D
	2.0	IPD/WB-IR-03/3D
	3.0	IPD/WB-IR-05/3D

Please see page 87 for the  Custom Ti Base System® accessories.



SCREWS

Coating	Screwdriver	Seat	ø4.0
-	1.20 Hex	Flat	IPD/WB-TR-00
TiN	1.20 Hex	Flat	IPD/WB-TR-00/TiN
	1.20 Hex	Flat	IPD/WB-TR-01

The IPD/WB-TR-01 model is an abutment extractor screw.



STRAIGHT MULTI-UNIT ABUTMENT

Angulation	Gingival Height	ø4.0
Straight	1.5	IPD/WB-MR-01
	2.5	IPD/WB-MR-02
	3.5	IPD/WB-MR-03
	4.0	IPD/WB-MR-04
	5.0	IPD/WB-MR-05

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	φ4.0
Analog	IPD/WB-AR-00
Digital Analog	IPD/WB-AR-00/3D
Bottom Screw	IPD/KA-TA-00
Side screw	IPD/KA-TA-01



OVERDENTURE RESTORATIONS

PSD SYSTEM ACCESSORIES

Indication	φ4.0
Analog	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00
Carrier	IPD/KA-CL-14



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	Height (mm)	φ4.0
	-	IPD/WB-PR-06
	-	IPD/WB-PR-07
Coping screw		IPD/WB-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Height (mm)	φ4.0
1.5	IPD/WB-DR-01
3.0	IPD/WB-DR-03
5.0	IPD/WB-DR-05



INSTRUMENTS AND ACCESSORIES

SCREWDRIVERS


Type	Lengths (mm)	φ4.0
1.20 Hex	18	IPD/PH-20-18
	25	IPD/PH-20-25
	32	IPD/PH-20-32
TPA	18	IPD/KA-CT-18
	25	IPD/KA-CT-25
	32	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	φ4.0
Clinical RA Shank Handle	16	IPD/MC-00-16
Laboratory RA Shank Handle	70	IPD/ML-00-70
RA Shank Ratchet	90	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90





Analogs/Digital Analogs



Compatibility NOBEL BIO CARE®

NOBEL REPLACE® SELECT



AC Serie NOBEL REPLACE® SELECT

[$\phi 3.5$ / $\phi 4.3$ / $\phi 5.0$]

ANALOG / DIGITAL ANALOG


Platform	Reference	Digital reference
$\phi 3.5$	IPD/AC-AN-00	IPD/AC-AN-00/3D
$\phi 4.3$	IPD/AC-AR-00	IPD/AC-AR-00/3D
$\phi 5.0$	IPD/AC-AW-00	IPD/AC-AW-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



Compatibility KLOCKNER®

SK2-NK2

ESSENTIAL CONE®



CA Serie SK2-NK2

[$\phi 4.2$]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
$\phi 4.2$	IPD/CA-AR-00	IPD/CA-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



CB Serie ESSENTIAL CONE®

[$\phi 4.5$]


ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
$\phi 4.5$	IPD/CB-AR-00	IPD/CB-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01





Compatibility KLOCKNER®



CC Serie KL™

[Ø3.5 / Ø4.1 / Ø5.1]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
Ø3.5	IPD/AA-AN-00	IPD/AA-AN-00/3D
Ø4.1	IPD/AA-AR-00	IPD/AA-AR-00/3D
Ø5.1	IPD/AA-AW-00	IPD/AA-AW-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



CD Serie VEGA®

[Ø3.5 / Ø4.0 / Ø4.5]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
Ø3.5	IPD/CD-AN-00	IPD/CD-AN-00/3D
Ø4.0 Ø4.5	IPD/CD-AR-00	IPD/CD-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



Compatibility STRAUMANN®



DD Serie SRA

[Ø3.5 / Ø4.6]


ANALOG / DIGITAL ANALOG

Platform	Reference	Digital Reference
Ø3.5	IPD/DD-AN-00	IPD/DD-AN-00/3D
Ø4.6	IPD/DD-AR-00	IPD/DD-AR-00/3D
Bottom screw	IPD/KA-TA-00	-
Side screw	IPD/KA-TA-01	-





Compatibility ZIMMER®



FB Serie SWISSPLUS

[ø4.8]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
ø4.8	IPD/DA-AR-00	IPD/DA-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



FC Serie EZTETIC®

[ø3.1]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
ø3.1	IPD/FC-AT-00	IPD/FC-AT-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



Compatibility BTI®



GA Serie EXTERNAL

[ø3.5 / ø4.1 / ø5.5]


ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
ø3.5 Tiny®	IPD/GA-AN-00	IPD/GA-AN-00/3D
ø4.1	IPD/AA-AR-00	IPD/AA-AR-00/3D
ø5.5	IPD/GA-AW-00	-
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01





Compatibility BTI®



GB Serie INTERNAL UNIVERSAL®

[Universal® $\phi 4.1$ / Wide $\phi 5.5$]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
$\phi 4.1$	IPD/GB-AR-00	IPD/GB-AR-00/3D
$\phi 5.5$	IPD/GB-AW-00	-
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



GC Serie MULTI-IM®

[$\phi 4.1$]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
$\phi 4.1$	IPD/GC-AR-00	IPD/GC-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



Compatibility MICRODENT®




HA Serie MICRODENT® SYSTEM

[$\phi 3.5$ / $\phi 4.2$ / $\phi 5.1$]

ANALOG

Platform	Reference
$\phi 3.5$	IPD/HA-AN-00
$\phi 4.2$	IPD/HA-AR-00
$\phi 5.1$	IPD/HA-AW-00





Compatibility MICRODENT®

MICRODENT® UNIVERSAL™



HB Serie MICRODENT® UNIVERSAL™

[ø3.5 / ø4.1 / ø5.1]

ANALOG / DIGITAL ANALOG


Platform	Reference	Digital reference
ø3.5	IPD/HA-AN-00	-
ø4.1	IPD/AA-AR-00	IPD/AA-AR-00/3D
ø5.1	IPD/BA-AW-00	IPD/BA-AW-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



Compatibility DENTSPLY®

XIVE® FRIADENT®

DENTSPLY®



IA Serie XIVE® FRIADENT®

[ø3.4 / ø3.8 / ø4.5]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
ø3.4	IPD/IA-AN-00	IPD/IA-AN-00/3D
ø3.8	IPD/IA-AR-00	IPD/IA-AR-00/3D
ø4.5	IPD/IA-AW-00	-
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



IB Serie ANKYLOS®

[C/X]


ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
C/X	IPD/IB-AR-00	IPD/IB-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01





Compatibility CAMLOG®



JA Serie CAMLOG®

[ø3.3 / ø3.8 / ø4.3]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
ø3.3	IPD/JA-AN-00	IPD/JA-AN-00/3D
ø3.8	IPD/JA-AR-00	IPD/JA-AR-00/3D
ø4.3	IPD/JA-AW-00	IPD/JA-AW-00/3D
ø5.0	IPD/JA-AX-00	IPD/JA-AX-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



JB Serie CONELOG®

[ø3.3 / ø3.8 / ø4.3]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
ø3.3	IPD/JB-AT-00	IPD/JB-AT-00/3D
ø3.8	IPD/JB-AN-00	IPD/JB-AN-00/3D
ø4.3	IPD/JB-AR-00	IPD/JB-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



Compatibility SWEDEN & MARTINA®




MA Serie OUTLINK®

[ø3.3 / ø4.1 / ø5.0]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
ø3.3	IPD/MA-AN-00	IPD/MA-AN-00/3D
ø4.1	IPD/AA-AR-00	IPD/AA-AR-00/3D
ø5.0	IPD/BA-AW-00	IPD/BA-AW-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01





Compatibility SWEDEN & MARTINA®

PREMIUM™ KHONO®



MB Serie PREMIUM™ KHONO®

[$\phi 3.3$ / $\phi 3.8$ / $\phi 4.25$ / $\phi 5.0$]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
$\phi 3.3$	IPD/MB-AT-00	IPD/MB-AT-00/3D
$\phi 3.8$	IPD/MB-AN-00	IPD/MB-AN-00/3D
$\phi 4.25$	IPD/MB-AR-00	IPD/MB-AR-00/3D
$\phi 5.0$	IPD/MB-AW-00	-
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



Compatibility PHIBO®

TSH®



NA Serie TSH®

[S2 / S3-4 / S5]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
S2	IPD/NA-AN-00	IPD/NA-AN-00/3D
S3-S4	IPD/AA-AR-00	IPD/AA-AR-00/3D
S5	IPD/HA-AW-00	IPD/NA-AW-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01





Compatibility BIOTECH[®] DENTAL

TAPERED INTERNAL



QB Serie KONTACT[®] [RP]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
RP	IPD/QB-AR-00	IPD/QB-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



Compatibility NEODENT[®]

HELIX[®] HE



RA Serie HELIX[®] HE [ø3.3 / ø4.1 / ø5.0]

ANALOG / DIGITAL ANALOG

Platform	Reference	Reference
ø3.3	IPD/RA-AN-00	IPD/RA-AN-00/3D
ø4.1	IPD/AA-AR-00	IPD/AA-AR-00/3D
ø5.0	IPD/RA-AW-00	IPD/RA-AW-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side screw	-	IPD/KA-TA-01





Compatibility BEGO® SEMADOS®

SC7RS



SB Serie SC/RS

[$\phi 3.25$ / $\phi 3.75$ / $\phi 4.1$ / $\phi 4.5$]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
$\phi 3.25$ $\phi 3.75$	IPD/SB-AN-00	IPD/SB-AN-00/3D
$\phi 4.1$	IPD/SB-AR-00	IPD/SB-AR-00/3D
$\phi 4.5$	IPD/SB-AW-00	IPD/SB-AW-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



Compatibility MIS®

V3®



TB Serie C1/V3®

[Standard]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
Standard	IPD/TB-AR-00	IPD/TB-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01





Compatibility DIO®



UB Serie UFII® [NP / RP]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
NP	IPD/UB-AN-00	IPD/UB-AN-00/3D
RP	IPD/UB-AR-00	IPD/UB-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



Compatibility DENTIUM®



XA Serie IMPLANTIUM® / SUPERLINE™ [RP]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
RP	IPD/XA-AR-00	IPD/XA-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01





Compatibility MEDENTIS®

ICX



YB Serie ICX

[RP]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
RP	IPD/YB-AR-00	IPD/YB-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



Compatibility GLOBAL D®

IN-KHONE®



1A Serie IN-KONE®

[RP]

ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
RP	IPD/1A-AR-00	IPD/1A-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01





Compatibility

ANTHOGYR®

AXIOME® BL



2A Serie AXIOME® BL

[RP]

ANALOG / DIGITAL ANALOG


Platform	Reference	Digital reference
RP	IPD/2A-AR-00	IPD/2A-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



Working Systems

-  Custom Ti Base System®
- MULTI-UNIT SYSTEM
- ACCESSORIES
- TECHNICAL SPECIFICATIONS

WORKFLOW DIAGRAM





CUTTING GUIDE

Engaging

Platform	Height (mm)	Reference
Narrow	3.5	IPD/KA-IN-35
	5.0	IPD/KA-IN-50
	6.5	IPD/KA-IN-65
Standard	3.5	IPD/KA-IR-35
	5.0	IPD/KA-IR-50
	6.5	IPD/KA-IR-65
Wide	3.5	IPD/KA-IW-35
	5.0	IPD/KA-IW-50
	6.5	IPD/KA-IW-65

Version with agulation window subject to availability.



NON-ENGAGING

Platform	Height (mm)	Reference
Narrow	3.5	IPD/KR-IN-35
	5.0	IPD/KR-IN-50
	6.5	IPD/KR-IN-65
Standard	3.5	IPD/KR-IR-35
	5.0	IPD/KR-IR-50
	6.5	IPD/KR-IR-65
Wide	3.5	IPD/KR-IW-35
	5.0	IPD/KR-IW-50
	6.5	IPD/KR-IW-65

Version with agulation window subject to availability.



FIXING SCREW

	Screwdriver	Reference
Every cutting guide	1.20 Hex	IPD/KA-TI-00



CUTTING GUIDE KITS

Platform	Heights (mm)		
Narrow	3.5 / 5.0 / 6.5	IPD/KT-IN-00	IPD/KT-IN-01
Standard	3.5 / 5.0 / 6.5	IPD/KT-IR-00	IPD/KT-IR-01
Wide	3.5 / 5.0 / 6.5	IPD/KT-IW-00	IPD/KT-IW-01

All kits include the locking screw IPD/KA-TI-00



SCAN BODY SCREWDRIVER

Scan Body height (mm)	Reference
10	IPD/PH-20-18
15	IPD/PH-20-25



Compatibility	Platform	Connection	GH (mm)	Reference	Cutting guides
Nobel Biocare® Branemark System®	Ø3.5	Engaging	1.5	IPD/AA-IN-02/3D	IPD/KT-IN-00
		Non-engaging		IPD/AA-IN-03/3D	IPD/KT-IN-01
	Ø4.1	Engaging	0.5	IPD/AA-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/AA-IR-02/3D	
			2.5	IPD/AA-IR-04/3D	
		Non-engaging	0.5	IPD/AA-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/AA-IR-03/3D	
			2.5	IPD/AA-IR-05/3D	
	Ø5.1	Engaging	1.3	IPD/AA-IW-02/3D	IPD/KT-IW-00
		Non-engaging		IPD/AA-IW-03/3D	-
Nobel Biocare® Multi-unit System	Ø4.8	Non-engaging	0.2	IPD/AB-IR-01/3D	-
			0.2	IPD/AB-IR-03/3D	
Nobel Biocare® Nobel Active®	Ø3.0	Engaging	0.8	IPD/AD-IT-00/3D	IPD/KT-IN-00
			1.5	IPD/AD-IT-02/3D	
			2.5	IPD/AD-IT-04/3D	
		Non-engaging	0.8	IPD/AD-IT-01/3D	IPD/KT-IN-01
			1.5	IPD/AD-IT-03/3D	
			2.5	IPD/AD-IT-05/3D	
	Ø3.5	Engaging	0.5	IPD/AD-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/AD-IN-02/3D	
			2.5	IPD/AD-IN-04/3D	
		Non-engaging	3.5	IPD/AD-IN-06/3D	IPD/KT-IN-01
			0.5	IPD/AD-IN-01/3D	
			1.5	IPD/AD-IN-03/3D	
	Ø4.3	Engaging	2.5	IPD/AD-IN-05/3D	IPD/KT-IR-00
			0.5	IPD/AD-IR-00/3D	
			1.5	IPD/AD-IR-02/3D	
		Non-engaging	2.5	IPD/AD-IR-04/3D	IPD/KT-IR-01
			0.5	IPD/AD-IR-01/3D	
			1.5	IPD/AD-IR-03/3D	
Biomet® 3i® Osseotite®	Ø3.4	Engaging	1.5	IPD/BA-IN-02/3D	IPD/KT-IN-00
		Non-engaging		IPD/BA-IN-03/3D	IPD/KT-IN-01
	Ø4.1	Engaging	0.5	IPD/AA-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/AA-IR-02/3D	
			2.5	IPD/AA-IR-04/3D	
	Ø5.0	Engaging	0.5	IPD/BA-IW-00/3D	IPD/KT-IW-00
			1.5	IPD/BA-IW-02/3D	
		Non-engaging	0.5	IPD/BA-IW-01/3D	IPD/KT-IW-01



Compatibility	Platform	Connection	GH (mm)	Reference	Cutting guides
Biomet® 3i® Certain®	Ø3.4	Engaging	0.5	IPD/BB-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/BB-IN-02/3D	
			2.5	IPD/BB-IN-04/3D	
		Non-engaging	0.5	IPD/BB-IN-01/3D	IPD/KT-IN-01
			1.5	IPD/BB-IN-03/3D	
			2.5	IPD/BB-IN-05/3D	
	Ø4.1	Engaging	0.5	IPD/BB-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/BB-IR-02/3D	
		Non-engaging	0.5	IPD/BB-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/BB-IR-03/3D	
	Ø5.0	Engaging	0.5	IPD/BB-IW-00/3D	IPD/KT-IW-00
			1.5	IPD/BB-IW-02/3D	
Non-engaging		0.5	IPD/BB-IW-01/3D	IPD/KT-IW-01	
		1.5	IPD/BB-IW-03/3D		
Straumann® Tissue Level	Ø4.8 RN	Engaging	0.4	IPD/DA-IR-00/3D	IPD/KT-IR-00
		Non-engaging		IPD/DA-IR-01/3D	IPD/KT-IR-01
	Ø6.5 WN	Engaging		IPD/DA-IW-02/3D	IPD/KT-IW-00
		Non-engaging		IPD/DA-IW-03/3D	IPD/KT-IW-01
		Straumann® Bone Level		Ø3.3 NC	Engaging
1.5	IPD/DB-IN-02/3D				
2.5	IPD/DB-IN-04/3D				
Non-engaging	3.5		IPD/DB-IN-06/3D		IPD/KT-IN-01
	0.5		IPD/DB-IN-01/3D		
	1.5		IPD/DB-IN-03/3D		
Ø4.1 RC	Engaging		2.5	IPD/DB-IN-05/3D	IPD/KT-IR-00
			3.5	IPD/DB-IN-07/3D	
			Non-engaging	0.5	IPD/DB-IR-00/3D
	1.5			IPD/DB-IR-02/3D	
	2.5			IPD/DB-IR-04/3D	
	Straumann® BLX®		RB/WB	Engaging	3.5
0.5		IPD/DB-IR-01/3D			
1.5		IPD/DB-IR-03/3D			
Non-engaging		2.5		IPD/DB-IR-05/3D	IPD/KT-IR-01
		3.5		IPD/DB-IR-07/3D	
Straumann® BLX®	RB/WB	Engaging	1.0	IPD/DC-IR-00/3D	IPD/KT-IR-00
			2.0	IPD/DC-IR-02/3D	
			3.0	IPD/DC-IR-04/3D	
		Non-engaging	1.0	IPD/DC-IR-01/3D	IPD/KT-IR-01
			2.0	IPD/DC-IR-03/3D	
			3.0	IPD/DC-IR-05/3D	

Compatibility	Platform	Connection	GH (mm)	Reference	Cutting guides
Astra® Osseospeed®	Ø3.0	Engaging	0.8	IPD/EA-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/EA-IN-02/3D	
			2.5	IPD/EA-IN-04/3D	
		Non-engaging	0.8	IPD/EA-IN-01/3D	IPD/KT-IN-01
			1.5	IPD/EA-IN-03/3D	
			2.5	IPD/EA-IN-05/3D	
	Ø3.5-4.0	Engaging	0.5	IPD/EA-IR-00/3D	IPD/KT-IN-00
			1.5	IPD/EA-IR-02/3D	
			2.5	IPD/EA-IR-04/3D	
		Non-engaging	3.5	IPD/EA-IR-06/3D	IPD/KT-IN-01
			0.5	IPD/EA-IR-01/3D	
			1.5	IPD/EA-IR-03/3D	
	Ø4.5-5.0	Engaging	2.5	IPD/EA-IR-05/3D	IPD/KT-IN-01
			3.5	IPD/EA-IR-07/3D	
			0.5	IPD/EA-IW-00/3D	IPD/KT-IR-00
		Non-engaging	1.5	IPD/EA-IW-02/3D	
			0.5	IPD/EA-IW-01/3D	
			1.5	IPD/EA-IW-03/3D	IPD/KT-IR-01
Astra® Evolution®	Ø3.6	Engaging	0.5	IPD/EB-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/EB-IN-02/3D	
			2.5	IPD/EB-IN-04/3D	
		Non-engaging	0.5	IPD/EB-IN-01/3D	IPD/KT-IN-01
			1.5	IPD/EB-IN-03/3D	
	Ø4.2	Engaging	2.5	IPD/EB-IN-05/3D	IPD/KT-IR-00
			0.5	IPD/EB-IR-00/3D	
			1.5	IPD/EB-IR-02/3D	
		Non-engaging	2.5	IPD/EB-IR-04/3D	IPD/KT-IR-01
			0.5	IPD/EB-IR-01/3D	
Zimmer® Screw Vent®	Ø3.5	Engaging	1.5	IPD/EB-IR-03/3D	IPD/KT-IR-01
			2.5	IPD/EB-IR-05/3D	
			0.5	IPD/FA-IN-00/3D	IPD/KT-IN-00
		Non-engaging	1.5	IPD/FA-IN-02/3D	
			2.5	IPD/FA-IN-04/3D	
	Ø4.5	Engaging	0.5	IPD/FA-IN-01/3D	IPD/KT-IN-01
			1.5	IPD/FA-IN-03/3D	
			2.5	IPD/FA-IN-05/3D	
		Non-engaging	0.5	IPD/FA-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/FA-IR-02/3D	
			2.5	IPD/FA-IR-04/3D	
	Ø5.7	Engaging	0.5	IPD/FA-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/FA-IR-03/3D	
		Non-engaging	2.5	IPD/FA-IR-05/3D	IPD/KT-IR-01
			0.5	IPD/FA-IW-00/3D	
	Ø5.7	Non-engaging	0.5	IPD/FA-IW-01/3D	IPD/KT-IW-00
			0.5	IPD/FA-IW-01/3D	


Compatibility	Platform	Connection	GH (mm)	Reference	Cutting guides
BioHorizons® Tapered Internal	Ø3.0	Engaging	0.5	IPD/LB-IT-00/3D	IPD/KT-IN-00
			1.5	IPD/LB-IT-02/3D	
			2.5	IPD/LB-IT-04/3D	
		Non-engaging	0.5	IPD/LB-IT-01/3D	IPD/KT-IN-01
			1.5	IPD/LB-IT-03/3D	
			2.5	IPD/LB-IT-05/3D	
	Ø3.5	Engaging	0.5	IPD/LB-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/LB-IN-02/3D	
			2.5	IPD/LB-IN-04/3D	
		Non-engaging	0.5	IPD/LB-IN-01/3D	IPD/KT-IN-01
			1.5	IPD/LB-IN-03/3D	
			2.5	IPD/LB-IN-05/3D	
	Ø4.5	Engaging	0.5	IPD/LB-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/LB-IR-02/3D	
			2.5	IPD/LB-IR-04/3D	
		Non-engaging	0.5	IPD/LB-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/LB-IR-03/3D	
			2.5	IPD/LB-IR-05/3D	
	Ø5.7	Engaging	0.5	IPD/LB-IW-00/3D	IPD/KT-IW-00
		Non-engaging	0.5	IPD/LB-IW-01/3D	IPD/KT-IW-01
Hiossen Implant® (Osstem) TSIII®	Ø3.5 Mini	Engaging	1.0	IPD/OB-IN-00/3D	IPD/KT-IN-00
			2.0	IPD/OB-IN-02/3D	
			3.0	IPD/OB-IN-04/3D	
		Non-engaging	1.0	IPD/OB-IN-01/3D	IPD/KT-IN-01
			2.0	IPD/OB-IN-03/3D	
			3.0	IPD/OB-IN-05/3D	
	Ø4.0 Regular	Engaging	1.0	IPD/OB-IR-00/3D	IPD/KT-IR-00
			2.0	IPD/OB-IR-02/3D	
			3.0	IPD/OB-IR-04/3D	
		Non-engaging	1.0	IPD/OB-IR-01/3D	IPD/KT-IR-01
			2.0	IPD/OB-IR-03/3D	
			3.0	IPD/OB-IR-05/3D	
Neodent® Grand Morse™ GM	GM	Engaging	1.0	IPD/RB-IR-00/3D	IPD/KT-IR-00
			2.0	IPD/RB-IR-02/3D	
			3.0	IPD/RB-IR-04/3D	
			4.0	IPD/RB-IR-06/3D	
		Non-engaging	1.0	IPD/RB-IR-01/3D	IPD/KT-IR-01
			2.0	IPD/RB-IR-03/3D	
			3.0	IPD/RB-IR-05/3D	
			4.0	IPD/RB-IR-07/3D	
		Engaging	0.5	IPD/RC-IR-00/3D	IPD/KT-IW-00
			0.5	IPD/RD-IN-00/3D	
Neodent® GM Abutment	RP	Engaging	0.5	IPD/RC-IR-00/3D	IPD/KT-IW-00
		Non-engaging	0.3	IPD/RD-IN-01/3D	IPD/KT-IN-01
Neodent® GM Micro Abutment	NP	Engaging	0.5	IPD/RD-IN-00/3D	-
		Non-engaging	0.3	IPD/RD-IN-01/3D	IPD/KT-IN-01

Compatibility	Platform	Connection	GH (mm)	Reference	Cutting guides
MIS® Seven®	Narrow Ø3.3	Engaging	0.5	IPD/TA-IT-00/3D	IPD/KT-IN-00
			1.5	IPD/TA-IT-02/3D	
			2.5	IPD/TA-IT-04/3D	
		Non-engaging	0.5	IPD/TA-IT-01/3D	IPD/KT-IN-01
			1.5	IPD/TA-IT-03/3D	
	Standard	Engaging	2.5	IPD/TA-IT-05/3D	IPD/KT-IN-00
			0.5	IPD/FA-IN-00/3D	
			1.5	IPD/FA-IN-02/3D	
		Non-engaging	2.5	IPD/FA-IN-04/3D	IPD/KT-IN-01
			0.5	IPD/FA-IN-01/3D	
	Wide	Engaging	1.5	IPD/FA-IN-03/3D	IPD/KT-IN-01
			2.5	IPD/FA-IN-05/3D	
			0.5	IPD/FA-IR-00/3D	IPD/KT-IR-00
		Non-engaging	1.5	IPD/FA-IR-02/3D	
			2.5	IPD/FA-IR-04/3D	
MIS® C1®	Standard	Engaging	0.5	IPD/FA-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/FA-IR-03/3D	
			2.5	IPD/FA-IR-05/3D	
			0.7	IPD/TB-IR-00/3D	IPD/KT-IR-00
		Non-engaging	1.5	IPD/TB-IR-02/3D	
			2.5	IPD/TB-IR-04/3D	
			3.5	IPD/TB-IR-06/3D	IPD/KT-IR-01
Megagen® AnyOne®	RP	Engaging	0.7	IPD/TB-IR-01/3D	
			1.5	IPD/TB-IR-03/3D	IPD/KT-IR-01
			2.5	IPD/TB-IR-05/3D	
		Non-engaging	3.5	IPD/TB-IR-07/3D	
			1.0	IPD/WA-IR-00/3D	IPD/KT-IR-00
Megagen® AnyRidge®	RP	Engaging	2.0	IPD/WA-IR-02/3D	
			3.0	IPD/WA-IR-04/3D	
			1.0	IPD/WA-IR-01/3D	IPD/KT-IR-01
		Non-engaging	2.0	IPD/WA-IR-03/3D	
			3.0	IPD/WA-IR-05/3D	
Megagen® AnyRidge®	RP	Engaging	1.0	IPD/WB-IR-00/3D	IPD/KT-IR-00
			2.0	IPD/WB-IR-02/3D	
			3.0	IPD/WB-IR-04/3D	
		Non-engaging	1.0	IPD/WB-IR-01/3D	IPD/KT-IR-01
			2.0	IPD/WB-IR-03/3D	
			3.0	IPD/WB-IR-05/3D	



WORKING SYSTEMS PSD System

WORKFLOW DIAGRAM



ANALOG

Description	Reference	Digital Reference
PSD analog	IPD/PS-AR-00	IPD/PS-AR-00/3D



TITANIUM INSERT

Description	Reference
Titanium insert	IPD/KA-CL-00
Titanium insert + 4 Retention inserts	IPD/KA-CL-02
Titanium insert + 4 Divergent Retention inserts	IPD/KA-CL-03



INSERTS

Description	Retention	Units	Reference
Combined inserts up to 20°	-	Pack 4	IPD/KA-CL-01
Red divergent inserts	1 lbs	Pack 4	IPD/KA-CL-04
Orange divergent inserts	2 lbs	Pack 4	IPD/KA-CL-05
Green divergent inserts	4 lbs	Pack 4	IPD/KA-CL-06
Inserts up to 20° in blue	1.5 lbs	Pack 8	IPD/KA-CL-10
Inserts up to 20° in grey	3 lbs	Pack 8	IPD/KA-CL-11
Inserts up to 20° colourless	5 lbs	Pack 8	IPD/KA-CL-12
Black processing inserts	-	Pack 8	IPD/KA-CL-13




ABUTMENT FOR BARS

Platform	Height (mm)	Metric	Units	Reference
PSD para barras	1	2 mm	2	IPD/PS-LR-00



PSD SCREWDRIVER

Description	Reference
	IPD/LL-PS-00



PSD TOOL


Description	Reference
PSD Screwdriver + insert extractor tool	IPD/LL-PS-01



WORKING SYSTEMS

Multi-unit System

WORKFLOW DIAGRAM



INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	ø4.8
	10	IPD/AB-SR-01

Captive screw included. 1.20 Hex Screwdriver.

SCAN TRANSFER

Indication	Height (mm)	Screw	ø4.8
	12.5	IPD/AB-TR-11/C	IPD/AB-SR-11

IPD/AB-SR-11 includes screw and a protective cap.

Ti-BASE Custom Ti-Base System®

Indication	Post Height (mm)	Gingival Height	ø4.8
	3.5	0.2	IPD/AB-IR-01/3D
	4.0	0.2	IPD/AB-IR-11/3D
	6.5	0.2	IPD/AB-IR-03/3D

Please see page 87 for the  Custom Ti Base System® accessories.

SCREWS

Coating	Screwdriver	Seat	Height	ø4.8
-	Unigrip	Flat	-	IPD/AB-TR-00
-			-	IPD/AB-TR-00/TIN
TiN	TPA	Flat	For 3.5 and 4.0 mm height ti-base	IPD/AB-TR-50
			For 6.5 mm height ti-base	IPD/AB-TR-51
-	TPA	Flat	-	IPD/EA-TR-51

IPD/EA-TR-51 is the dedicated screw for MUA level restorations on Zr-PMMA.

ANALOGS

Types	ø4.8
Analog	IPD/AB-AR-00
Digital Analog	IPD/AB-AR-00/3D
Bottom Screw	IPD/KA-TA-00
Side Screw	IPD/KA-TA-01

TEMPORARY ABUTMENT / OPEN TRAY COPING

Indication	ø4.8
	IPD/AB-PR-01
	IPD/AB-PR-00
	IPD/AB-PR-06
Coping screw	IPD/AB-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.
IPD/AB-PR-06 corresponds to the Classic design and cannot be used as a coping.



WORKING SYSTEMS

Multi-unit System

OPEN AND CLOSED TRAY COPING

	ø4.8
Open tray coping	IPD/AB-CR-00
Closed tray coping	IPD/AB-CR-01



HEALING ABUTMENT

Material	ø4.8
Plastic POM	IPD/AB-DR-00
Titanium	IPD/AB-DR-01

IPD/AB-TR-00 screw is not included.



PSD FOR MULTI-UNIT

Units	Height (mm)	ø4.8
2	3.35	IPD/AB-LR-02



PSD System Accessories

Indication	ø4.8
Analog	IPD/PS-AR-00
Digital Analog	IPD/PS-AR-00/3D
Bottom screw	IPD/KA-TA-00
Carrier	IPD/KA-CL-14



ANGLED MULTI-UNIT CARRIER

Indication	ø4.8
	IPD/LL-TM-00



SCREWDRIVERS

Type	Lengths (mm)	ø4.8
Unigrip	18	IPD/PU-00-18
	25	IPD/PU-00-25
	32	IPD/PU-00-32
2.0 Hex	Straight Multi-Unit abutment	IPD/LL-AB-00



WORKING SYSTEMS

Multi-unit System

Compatibility	Platform	Angulation	Height (mm)	Reference
Nobel Biocare® Branemark System®	ø3.5	Straight	1.0	IPD/AA-MR-01
			2.0	IPD/AA-MR-02
			3.0	IPD/AA-MR-03
	ø4.1	Straight	1.0	IPD/AB-MR-01
			2.0	IPD/AB-MR-02
			3.0	IPD/AB-MR-03
			4.0	IPD/AB-MR-04
			5.0	IPD/AB-MR-05
	ø3.5	Straight	1.5	IPD/AD-MN-01
			2.5	IPD/AD-MN-02
			3.5	IPD/AD-MN-03
			4.5	IPD/AD-MN-04
		17°	2.5	IPD/AD-MN-12
			3.5	IPD/AD-MN-13
		30°	3.5	IPD/AD-MN-33
			4.5	IPD/AD-MN-34
			4.5	IPD/AD-MN-34
Nobel Biocare® Nobel Active®	ø4.3	Straight	1.5	IPD/AD-MR-01
			2.5	IPD/AD-MR-02
			3.5	IPD/AD-MR-03
			4.5	IPD/AD-MR-04
		17°	2.5	IPD/AD-MR-12
			3.5	IPD/AD-MR-13
	ø4.3	30°	3.5	IPD/AD-MR-33
			4.5	IPD/AD-MR-34
			4.5	IPD/AD-MR-34
		Straight	1.0	IPD/AB-MR-01
			2.0	IPD/AB-MR-02
Biomet® 3i® Osseotite®	ø4.1	Straight	3.0	IPD/AB-MR-03
			4.0	IPD/AB-MR-04
			5.0	IPD/AB-MR-05
			5.0	IPD/AB-MR-05
			5.0	IPD/AB-MR-05
Straumann® Tissue Level	ø4.8 RN	Straight	0.9	IPD/DA-MR-01
			1.5	IPD/DC-MR-01
Straumann® BLX®	RB/WB	Straight	2.5	IPD/DC-MR-02
			3.5	IPD/DC-MR-03
		17°	3.5	IPD/DC-MR-13
			4.5	IPD/DC-MR-34
		30°	4.5	IPD/DC-MR-34

WORKING SYSTEMS

Multi-unit System

Compatibility	Platform	Angulation	Height (mm)	Reference
Straumann® Bone Level	ø3.3 NC	Straight	1.5	IPD/DB-MN-01
			2.5	IPD/DB-MN-02
			3.5	IPD/DB-MN-03
			4.5	IPD/DB-MN-04
		17°	2.5	IPD/DB-MN-12
			3.5	IPD/DB-MN-13
		30°	3.5	IPD/DB-MN-33
			4.5	IPD/DB-MN-34
	ø4.1 RC	Straight	1.5	IPD/DB-MR-01
			2.5	IPD/DB-MR-02
			3.5	IPD/DB-MR-03
			4.5	IPD/DB-MR-04
		17°	2.5	IPD/DB-MR-12
			3.5	IPD/DB-MR-13
		30°	3.5	IPD/DB-MR-33
			4.5	IPD/DB-MR-34
Astra® Osseospeed™	ø3.5 / ø4.0	Straight	1.5	IPD/EA-MR-01
			2.5	IPD/EA-MR-02
			3.5	IPD/EA-MR-03
			4.5	IPD/EA-MR-04
		17°	2.5	IPD/EA-MR-12
			3.5	IPD/EA-MR-13
		30°	3.5	IPD/EA-MR-33
			4.5	IPD/EA-MR-34
	ø4.5 / ø5.0	Straight	1.5	IPD/EA-MW-01
			2.5	IPD/EA-MW-02
			3.5	IPD/EA-MW-03
			4.5	IPD/EA-MW-04
		17°	2.5	IPD/EA-MW-12
			3.5	IPD/EA-MW-33
Astra® Evolution®	ø3.6	Straight	1.5	IPD/EB-MN-01
			2.5	IPD/EB-MN-02
			3.5	IPD/EB-MN-03
		17°	2.5	IPD/EB-MN-12
			3.5	IPD/EB-MN-33
			3.5	IPD/EB-MN-33
	ø4.2	Straight	1.5	IPD/EB-MR-01
			2.5	IPD/EB-MR-02
			3.5	IPD/EB-MR-03
		17°	2.5	IPD/EB-MR-12
			3.5	IPD/EB-MR-33
			3.5	IPD/EB-MR-33
Zimmer® Screw Vent®	ø3.5	Straight	0.9	IPD/FA-MN-00
			1.5	IPD/FA-MN-01
			2.5	IPD/FA-MN-02
			3.5	IPD/FA-MN-03
			4.5	IPD/FA-MN-04
			4.5	IPD/FA-MN-04

WORKING SYSTEMS

Multi-unit System

Compatibility	Platform	Angulation	Height (mm)	Reference
Zimmer® Screw Vent®	Ø4.5	Straight	1.5	IPD/FA-MR-01
			2.5	IPD/FA-MR-02
			3.5	IPD/FA-MR-03
			4.5	IPD/FA-MR-04
BioHorizons® Tapered Internal	Ø3.5	17°	2.5	IPD/LB-MN-12
			3.5	IPD/LB-MN-13
		30°	3.5	IPD/LB-MN-33
			4.5	IPD/LB-MN-34
Hiossen Implant® TSIII	Ø3.5	Straight	1.5	IPD/OB-MN-01
			2.5	IPD/OB-MN-02
			3.5	IPD/OB-MN-03
			4.5	IPD/OB-MN-04
		17°	5.5	IPD/OB-MN-05
			2.5	IPD/OB-MN-12
			3.5	IPD/OB-MN-13
			3.5	IPD/OB-MN-33
		30°	4.5	IPD/OB-MN-34
			1.5	IPD/OB-MR-01
			2.5	IPD/OB-MR-02
			3.5	IPD/OB-MR-03
Neodent® Grand Morse™ GM	GM	Straight	4.5	IPD/OB-MR-04
			5.5	IPD/OB-MR-05
			2.5	IPD/OB-MR-12
			3.5	IPD/OB-MR-13
		17°	3.5	IPD/OB-MR-33
			4.5	IPD/OB-MR-34
		30°	1.5	IPD/RB-MR-01
			2.5	IPD/RB-MR-02
MIS® Seven®	Ø3.3	Straight	3.5	IPD/RB-MR-03
			2.5	IPD/RB-MR-12
			3.5	IPD/RB-MR-33
			4.5	IPD/RB-MR-34
MIS® C1®	Standard	Straight	2.0	IPD/TA-MT-02
			3.0	IPD/TA-MT-03
			4.0	IPD/TA-MT-04
			1.0	IPD/TB-MR-01
Megagen® AnyOne®	RP	Straight	2.0	IPD/TB-MR-02
			3.0	IPD/TB-MR-03
			4.0	IPD/TB-MR-04
			1.5	IPD/WA-MR-01
Megagen® AnyRidge®	RP	Straight	2.5	IPD/WA-MR-02
			3.5	IPD/WA-MR-03
			1.5	IPD/WB-MR-01
			2.5	IPD/WB-MR-02
			3.5	IPD/WB-MR-03

Accessories

1.00 HEX SCREWDRIVER

Lengths (mm)	Reference
18	IPD/PH-00-18
25	IPD/PH-00-25
32	IPD/PH-00-32



1.18 HEX SCREWDRIVER

Lengths (mm)	Reference
18	IPD/PH-18-18
25	IPD/PH-18-25
32	IPD/PH-18-32



1.20 HEX SCREWDRIVER

Lengths (mm)	Reference
18	IPD/PH-20-18
25	IPD/PH-20-25
32	IPD/PH-20-32



1.25 HEX SCREWDRIVER

Lengths (mm)	Reference
18	IPD/PH-25-18
25	IPD/PH-25-25
32	IPD/PH-25-32



1.27 HEX SCREWDRIVER

Lengths (mm)	Reference
18	IPD/PH-27-18
25	IPD/PH-27-25
32	IPD/PH-27-32



1.40 HEX SCREWDRIVER

Lengths (mm)	Reference
18	IPD/PH-40-18
25	IPD/PH-40-25
32	IPD/PH-40-32



1.70 HEX SCREWDRIVER

Lengths (mm)	Reference
18	IPD/PH-70-18
25	IPD/PH-70-25
32	IPD/PH-70-32



1.25 SQUARE SCREWDRIVER

Lengths (mm)	Reference
18	IPD/PC-25-18
25	IPD/PC-25-25
32	IPD/PC-25-32



STAR SCREWDRIVER (Klockner Essential® / Vega®)

Lengths (mm)	Reference
18	IPD/PS-00-18
25	IPD/PS-00-25
32	IPD/PS-00-32



UNIGRIP SCREWDRIVER

Lengths (mm)	Reference
18	IPD/PU-00-18
25	IPD/PU-00-25
32	IPD/PU-00-32



NEO GM SCREWDRIVER

Lengths (mm)	Reference
18	IPD/PU-01-18
25	IPD/PU-01-25
32	IPD/PU-01-32



FLAT SCREWDRIVER

Lengths (mm)	Reference
18	IPD/PP-00-18
25	IPD/PP-00-25
32	IPD/PP-00-32



TORX SCREWDRIVER

Lengths (mm)	Reference
18	IPD/PT-00-18
25	IPD/PT-00-25
32	IPD/PT-00-32




MULTI-UNIT SCREWDRIVER

Lengths (mm)	Reference
20	IPD/LL-AB-00



PSD SCREWDRIVER

Lengths (mm)	Reference
25	IPD/LL-PS-00



Accessories

PSD AND INSERT EXTRACTOR TOOL

Reference
IPD/LL-PS-01



TPA SCREWDRIVER

Lengths (mm)	Reference
18	IPD/KA-CT-18
25	IPD/KA-CT-25
32	IPD/KA-CT-32



CLINICAL RA SHANK HANDLE

Lengths (mm)	Reference
16	IPD/MC-00-16



LABORATORY RA SHANK HANDLE

Lengths (mm)	Reference
70	IPD/ML-00-70



RA SHANK RATCHET

Lengths (mm)	Reference
90	IPD/DI-00-90

Ratchet adapter included.



RATCHET ADAPTER

Lengths (mm)	Reference
Ratchet adapter	IPD/DI-KA-90



Accessories kits

Description	Reference	
	IPD/KT-DI-18	IPD/KT-DI-25
	Lengths	
1.20 Hex Screwdriver	18 mm	25 mm
1.25 Hex Screwdriver	18 mm	25 mm
1.27 Hex Screwdriver	18 mm	25 mm
1.70 Hex Screwdriver	18 mm	25 mm
1.25 Square screwdriver	18 mm	25 mm
Unigrip screwdriver	18 mm	25 mm
Flat screwdriver	18 mm	25 mm
TPA screwdriver	18 mm	25 mm
Torx screwdriver	18 mm	25 mm
Multi-unit screwdriver	20 mm	20 mm
RA Shank Ratchet	90 mm	90 mm



Description	Reference	
	IPD/KT-MC-18	IPD/KT-MC-25
	Lengths	
1.20 Hex Screwdriver	18 mm	25 mm
1.25 Hex Screwdriver	18 mm	25 mm
1.27 Hex Screwdriver	18 mm	25 mm
1.70 Hex Screwdriver	18 mm	25 mm
1.25 Square screwdriver	18 mm	25 mm
Unigrip screwdriver	18 mm	25 mm
Flat screwdriver	18 mm	25 mm
TPA screwdriver	18 mm	25 mm
Torx screwdriver	18 mm	25 mm
Multi-unit screwdriver	20 mm	20 mm
Clinical RA Shank Handle	16 mm	16 mm



Description	Reference	
	IPD/KT-ML-32	
	Lengths	
1.20 Hex Screwdriver	32 mm	
1.25 Hex Screwdriver	32 mm	
1.27 Hex Screwdriver	32 mm	
1.70 Hex Screwdriver	32 mm	
1.25 Square screwdriver	32 mm	
Unigrip screwdriver	32 mm	
Flat screwdriver	32 mm	
TPA screwdriver	32 mm	
Torx screwdriver	32 mm	
Laboratory RA Shank Handle	70 mm	



Screws technical specifications

DESCRIPTION				CHARACTERISTIC							
REFERENCE	TiN Reference	COMPATIBILITY	PLATFORM	SEAT	HEAD DIAMETER	METRIC	HEAD (L1)	BODY (L2)	TOTAL (LT)	SCREW-DRIVER	TORQUE (NCM)
IPD/AA-TN-00	IPD/AA-TN-00/TIN	Nobel Biocare® Branemark System®	Ø3.5	70°	2.25	M1.6x0.35	2.15	5.15	7.30	Unigrip	35
		Nobel Biocare® Nobel Active®	Ø3.5								
IPD/AA-TR-00	IPD/AA-TR-00/TIN	Nobel Biocare® Branemark System®	Ø4.1	75°	2.4	M2x0.4	1.5	5.85	7.35	Unigrip	35
		Nobel Biocare® Nobel Active®	Ø4.3								
IPD/AA-TR-02	IPD/AA-TR-02/TIN	Nobel Biocare® Branemark System®	Ø4.1	Flat	2.40	M2x0.4	1.45	5.9	7.35	Unigrip	35
		Nobel Biocare® Nobel Active®	Ø4.3								
IPD/AA-TW-00	IPD/AA-TW-00/TIN	Nobel Biocare® Branemark System®	Ø5.1	72°	3.40	M2.5x0.45	1.37	5.33	6.7	Unigrip	35
IPD/AB-TR-00	IPD/AB-TR-00/TIN	Multi-unit	Ø4.8	Flat	2.00	M1.4x0.3	1.35	2.1	3.45	Unigrip	15
IPD/AD-TT-00	IPD/AD-TT-00/TIN	Nobel Biocare® Nobel Active®	Ø3.0	70°	1.82	M1.4x0.3	1.86	4.84	6.7	Unigrip	15
IPD/AA-TR-01	IPD/AA-TR-01/TIN	Biomet® 3i® Osseotite® (Plat. R)	Ø3.4 / Ø4.1 / Ø5.0	75°	2.5	M2x0.4	1.67	4.6	6.27	1.20 Hex	20
IPD/BA-TR-00	IPD/BA-TR-00/TIN	Biomet® 3i® Osseotite® (Plat. R)	Ø4.1	75°	2.35	M2x0.4	1.79	5.81	7.6	1.20 Hex	20
IPD/BA-TR-01	IPD/BA-TR-01/TIN	Biomet® 3i® Osseotite® (Plat. R)	Ø4.1	75°	2.35	M2x0.4	1.79	5.81	7.6	1.25 Hex	25
IPD/BA-TR-02	IPD/BA-TR-02/TIN	Biomet® 3i® Osseotite®	Ø3.4 / Ø4.1 / Ø5.0	Flat	2.4	M2x0.4	1.65	6.15	7.8	1.25 Square	35
IPD/BB-TR-01	IPD/BB-TR-01/TIN	Biomet® 3i® Certain®	Ø3.4 / Ø4.1 / Ø5.0	Flat	2.0	M1.6x0.35	1.9	6.5	8.4	1.20 Hex	20
IPD/BB-TR-02	-	Biomet® 3i® Certain®	Ø4.1 / Ø5.0	140°	2.30	M1.6x0.35	2.05	6.35	8.4	1.20 Hex	20
IPD/DA-TR-00	IPD/DA-TR-00/TIN	Straumann® Tissue Level (synOcta® level)	Ø4.8 RN	90°	2.4	M2x0.4	1.2	3.3	4.5	Torx	35
IPD/DA-TR-01	IPD/DA-TR-01/TIN	Straumann® Tissue Level (Implant level)	Ø4.8 RN / Ø6.5 WN	30°	2.4	M2x0.4	3.03	3.67	6.7	Torx	35
IPD/DA-TR-02	-	Straumann® Tissue Level	Ø4.8 RN	Flat	2.52	M2x0.4	1.45	5.9	7.35	Torx	35
IPD/DA-TR-03	IPD/DA-TR-03/TIN	Straumann® Tissue Level	Ø4.8 RN	30°	2.57	M2x0.4	3.05	3.65	6.7	Torx	35
IPD/DB-TR-00	IPD/DB-TR-00/TIN	Straumann® Bone Level	Ø3.3 NC / Ø4.1 RC	30°	2.1	M1.6x0.35	2.31	5.54	7.85	Torx	35
IPD/DC-TR-00	IPD/DC-TR-00/TIN	Straumann® BLX®	RB/WB	90°	2.1	M1.6x0.35	1.5	4.6	6.1	Torx	30
IPD/EA-TN-00	IPD/EA-TN-00/TIN	Astra® Osseospeed™	Ø3.0	120°	2.1	M1.4x0.3	2	5.5	7.5	1.27 Hex	15
IPD/EA-TR-00	IPD/EA-TR-00/TIN	Astra® Osseospeed™	Ø3.5 - Ø4.0	Flat	2.35	M1.6x0.35	1.95	6.3	8.25	1.27 Hex	20
IPD/EA-TR-01	IPD/EA-TR-01/TIN	Astra® Osseospeed™	Ø3.5-Ø4.0 / Ø4.5-Ø5.0	Flat	2.25	M1.4x0.3	1.85	2.95	4.8	1.27 Hex	15
IPD/EA-TW-00	IPD/EA-TW-00/TIN	Astra® Osseospeed™	Ø4.5 - Ø5.0	Flat	2.35	M2x0.4	1.9	8.4	10.3	1.27 Hex	25
IPD/EB-TN-00	IPD/EB-TN-00/TIN	Astra® Evolution®	Ø3.6	60°	2.15	M1.6x0.35	2.12	5.33	7.45	1.27 Hex	20
IPD/EB-TR-00	IPD/EB-TR-00/TIN	Astra® Evolution®	Ø4.2	100°	2.35	M1.8x0.35	1.92	5.70	7.62	1.27 Hex	20
IPD/FA-TR-00	IPD/FA-TR-00/TIN	Zimmer® Screw Vent®	Ø3.5/ Ø4.5 / Ø5.7	140°	2.1	M1.8x0.35	2.54	5.66	8.2	1.25 Hex	30
IPD/LB-TR-00	IPD/LB-TR-00/TIN	BioHorizons® Tapered Internal	Ø3 / Ø3.5 / Ø4.5 / Ø5.7	90°	2.35	M1.8x0.35	2.42	4.83	7.25	1.27 Hex	35
IPD/OB-TN-00	IPD/OB-TN-00/TIN	Hiossen Implant® (Osstem) TSIII	Ø3.5	90°	2.20	M1.6x0.35	2.41	7.79	10.2	1.20 Hex	20
IPD/OB-TR-00	IPD/OB-TR-00/TIN	Hiossen Implant® (Osstem) TSIII	Ø4.0	90°	2.30	M2x0.4	2.23	6.12	8.35	1.20 Hex	30
IPD/RB-TR-00	IPD/RB-TR-00/TIN	Neodent® Grand Morse™ GM	GM	60°	1.86	M1.6x0.35	1.84	7.93	9.77	Neo GM	20
IPD/RC-TR-00	IPD/RC-TR-00/TIN	Neodent® GM abutment	RP	90°	2.35	M2x0.25	1.53	2.25	3.78	Neo GM	10
IPD/RD-TN-00	IPD/RD-TN-00/TIN	Neodent® GM micro abutment	Ø3.5	Flat	1.90	M1.4x0.3	1.72	1.99	3.71	Neo GM	10
IPD/TA-TT-00	IPD/TA-TT-00/TIN	MIS® Seven®	Narrow 3.3	120°	2.1	M1.6x0.35	2.58	5.1	7.68	1.25 Hex	30
IPD/TB-TR-00	IPD/TB-TR-00/TIN	MIS® C1®	Standard	120°	2.2	UNF 1-72	2.44	5.18	7.62	1.25 Hex	30
IPD/WA-TR-00	IPD/WA-TR-00/TIN	Megagen® AnyOne®	RP	120°	2.3	M2x0.4	1.9	8.0	9.9	1.20 Hex	35
IPD/WB-TR-00	IPD/WB-TR-00/TIN	Megagen® AnyRidge®	RP	Flat	2.1	M1.8x0.35	2	6	8	1.20 Hex	25
IPD/VA-TR-00	-			Flat	2.4	M2x0.4	2.2	5.8	8	1.20 Hex	25
IPD/XA-TR-02	-	Nobel Biocare® Branemark System®	Ø4.1	Flat	2.52	M2x0.4	1.45	5.9	7.35	1.20 Hex	30
		Nobel Biocare® Nobel Active®	Ø4.3								

TPA screw technical specifications

DESCRIPTION			CHARCATERISTIC							
REFERENCE	COMPATIBILITY	PLATFORM	SEAT	HEAD DIAMETER	METRIC	HEAD (L1)	BODY (L2)	TOTAL (LT)	SCREW-DRIVER	TORQUE (NCM)
IPD/AA-TN-50	Nobel Biocare® Branemark System®	Ø3.5	70°	2.3	M1.6x0.35	2.08	4.72	6.8	TPA	20
	Nobel Biocare® Nobel Active®	Ø3.5								
IPD/AA-TR-50	Nobel Biocare® Branemark System®	Ø4.1	75°	2.3	M2x0.4	1.93	3.87	5.8	TPA	25
	Biomet® 3i® Osseotite®	Ø4.1	75°	2.3	M2x0.4	1.93	3.87	5.8	TPA	25
IPD/AA-TR-52	Nobel Biocare® Branemark System®	Ø4.1	Flat	2.3	M2x0.4	1.4	4.4	5.8	TPA	25
	Biomet® 3i® Osseotite®	Ø3.25 / Ø4.1 / Ø5.1								
IPD/AA-TW-50	Nobel Biocare® Branemark System®	Ø5.1	72°	2.85	M2.5x0.45	1.61	4.8	6.41	TPA	25
IPD/AB-TR-50	Multi-unit / Multi-im®	Ø4.8	Flat	2.3	M1.4x0.3	1.5	2.1	3.6	TPA	15
IPD/AB-TR-51	Multi-unit	Ø4.8	Flat	2.25	M1.4X0.3	3.25	2.10	5.35	TPA	15
IPD/AC-TR-50	Astra® OsseoSpeed®	Ø4.5-Ø5.0	Flat	2.3	M2x0.4	1.4	7.4	8.8	TPA	25
IPD/AC-TR-51	Astra® OsseoSpeed®	Ø4.3	Flat	2.3	M2x0.4	1.3	6.8	8.1	TPA	25
IPD/AD-TR-50	Nobel Biocare® Nobel Active®	Ø4.3	75°	2.3	M2x0.4	2.0	5.0	7	TPA	25
IPD/AD-TT-50	Nobel Biocare® Nobel Active®	Ø3.0	70°	2.3	M1.4x0.3	2.0	4.65	6.65	TPA	15
IPD/BB-TR-51	Biomet® 3i® Certain®	Ø3.4 / Ø4.1 / Ø5.0	Flat	2.25	M1.6x0.35	1.35	6.25	7.6	TPA	20
IPD/DA-TR-50	Straumann® Tissue Level (SynOcta®)	Ø4.8 RN / Ø6.5 WN	90°	2.3	M2x0.4	1.65	3.32	4.97	TPA	25
IPD/DA-TR-51	Straumann® Tissue Level (Implant level)	Ø4.8 RN / Ø6.5 WN	30°	2.3	M2x0.4	2.65	3.65	6.3	TPA	25
IPD/DB-TR-50	Straumann® Bone Level	Ø3.3 NC / Ø4.1 RC	30°	2.3	M1.6x0.35	2.15	5.55	7.7	TPA	20
IPD/EA-TN-50	Astra® OsseoSpeed®	Ø3.0	120°	2.25	M1.4x0.35	1.65	5.45	7.1	TPA	15
IPD/EA-TR-50	Astra® OsseoSpeed®	Ø3.5-Ø4.0	Flat	2.25	M1.6x0.35	1.35	6.05	7.4	TPA	20
IPD/EA-TR-51	Astra® OsseoSpeed®	All	Flat	2.3	M1.4x0.3	1.50	2.90	4.4	TPA	15
IPD/EB-TN-50	Astra® Evolution®	Ø3.6	60°	2.3	M1.6x0.35	1.99	5.11	7.1	TPA	20
IPD/EB-TR-50	Astra® Evolution®	Ø4.2	100°	2.3	M1.8x0.35	2.22	5.38	7.6	TPA	25
IPD/FA-TR-50	Zimmer® Screw Vent®	Ø3.5 / Ø4.5 / Ø5.7	140°	2.25	M1.8x0.35	1.73	4.97	6.7	TPA	25
IPD/LB-TR-50	BioHorizons® Tapered Internal	Ø3.0 / Ø3.5 / Ø4.5 / Ø5.7	90°	2.3	M1.8x0.35	1.89	4.61	6.5	TPA	25
IPD/OB-TN-50	Hiossen Implant® (Osstem) TSIII	Ø3.5	60°	2.3	M1.55x0.35	1.98	7.52	9.50	TPA	20
IPD/OB-TR-50	Hiossen Implant® (Osstem) TSIII	Ø4.0	90°	2.3	M2x0.4	1.36	6.14	7.50	TPA	25
IPD/RB-TR-50	Neodent® Grand Morse™ GM	GM	60°	2.3	M1.6x0.35	1.98	7.52	9.50	TPA	20
IPD/RC-TR-50	Neodent® GM abutment	RP	90°	2.35	M2x0.25	1.60	2.23	3.83	TPA	10
IPD/RD-TN-50	Neodent® GM micro abutment	Ø3.5	Flat	2.25	M1.4x0.3	1.40	1.99	3.99	TPA	10
IPD/TA-TT-50	MIS® Seven®	Narrow 3.3	120°	2.3	M1.6x0.35	1.72	4.88	6.6	TPA	25
IPD/TB-TR-50	MIS® C1®	Standard	120°	2.3	UNF 1-72	2.35	4.85	7.2	TPA	25
IPD/WA-TR-50	Megagen® AnyOne®	RP	120°	2.3	M2x0.4	1.8	7.95	9.75	TPA	25



Now you have seen our
complete range of products,
be part of *ipd*.

www.ipd2004.com





Notes



US Office

1330 SW 22nd Street, Suite 200
Miami, FL 33145

Phone: +1 305 392 1427
burbina@ipd2004.com

Spain Office

C/ Rosa dels Vents, 9-15
08338 Premià de Dalt (Barcelona)

Phone: (+34) 93 278 84 91
info@ipd2004.com
www.ipd2004.com

